





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)			HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION		
				TOTAL	THREAD	A B C							
AB I2 STANDARD PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
AB I22 STANDARD PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
AB I5 STANDARD PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
AB I55 STANDARD PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
AB I10 STANDARD PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ACE INFINITY EXTERNAL HEX NARROW PLATFORM	TPDH16L6	1,6	20	5,9	3	4,4		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ACE INFINITY EXTERNAL HEX REGULAR PLATFORM	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ACE INFINITY EXTERNAL HEX WIDE PLATFORM	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ACE INFINITY TRI-CAM NARROW PLATFORM	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
ACE INFINITY TRI-CAM REGULAR PLATFORM	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ACE INFINITY TRI-CAM WIDE PLATFORM	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ACE INFINITY INTERNAL HEX 3,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ACE INFINITY INTERNAL HEX 4,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN SWELL INTERNAL HEX Ø3,3	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN SWELL INTERNAL HEX Ø3,75	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN SWELL INTERNAL HEX Ø4,2	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN SWELL INTERNAL HEX Ø5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN SWELL INTERNAL HEX Ø6	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG-OS Ø3,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG-OS Ø3,75	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG-OS Ø4,2	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG-OS Ø5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG-OS Ø6	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG-S Ø3,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG-S Ø3,75	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG-S Ø4,2	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG-S Ø5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG-S Ø6	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG CLOSFIT Ø3,5	TPDH16L73	1,6	20	7,3	2,2	4,87	5,56	1,74	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG CLOSFIT Ø4,3	TPDH2L75	2	25	7,5	2,75	5,93	6,18	1,32	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
ADIN TOUAREG CLOSFIT Ø5	TPDH2L75	2	25	7,5	2,75	5,93	6,18	1,32	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
ALPHA BIO INTERNAL HEXAGON NARROW PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ALPHA BIO INTERNAL HEXAGON REGULAR PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ALPHA BIO INTERNAL HEXAGON WIDE PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ALPHA BIO EXTERNAL HEXAGON NARROW PLATFORM	TPDH16L6	1,6	20	5,9	3	4,4		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ALPHA BIO EXTERNAL HEXAGON REGULAR PLATFORM	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ALPHA BIO EXTERNAL HEXAGON WIDE PLATFORM	TPDH25L67	2,5	25	6,7	4,6	5,1		1,6	2,85	Straight		Semiespherical	1,70 Hexalobular
ANCLADEN ANCLALOCK Ø3,75	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ANCLADEN ANCLALOCK Ø4,25	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ANCLADEN ANCLALOCK Ø5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ANTHOGRYR AXIOM REG Ø3,4	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ANTHOGRYR AXIOM REG Ø4	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ANTHOGRYR AXIOM REG Ø4,6	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ANTHOGRYR AXIOM REG Ø5,2	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ANTHOGRYR AXIOM PX Ø3,4	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ANTHOGRYR AXIOM PX Ø4	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ANTHOGRYR AXIOM PX Ø4,6	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ANTHOGRYR AXIOM PX Ø5,2	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
ANTHOGRYR OSSFIT PLATFORM Ø4,8	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)			HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION		
				TOTAL	THREAD	A						B	C
ANTHOGRYR OSSFIT PLATFORM Ø6,5	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular	
ANTHOGRYR ANTHOFIT HE EXTERNAL HEXAGON Ø4,1	TPDH2L74	2	25	7,4	3,3	6			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
ANTHOGRYR ANTHOFIT HE EXTERNAL HEXAGON Ø5	TPDH2L74	2	25	7,4	3,3	6			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
ANTHOGRYR ANTHOFIT OI INTERNAL OCTAGON Ø3,5	TPDH2L74	2	25	7,4	3,3	6			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
ANTHOGRYR ANTHOFIT OI INTERNAL OCTAGON Ø4	TPDH2L74	2	25	7,4	3,3	6			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
ANTHOGRYR ANTHOFIT OI INTERNAL OCTAGON Ø5	TPDH2L74	2	25	7,4	3,3	6			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
ARDS SMART HEXÁGONO INTERNO Ø3,75	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS SMART HEXÁGONO INTERNO Ø4,2	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS SMART HEXÁGONO INTERNO Ø4,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS CLASSIC HEXÁGONO INTERNO Ø3,3	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS CLASSIC HEXÁGONO INTERNO Ø3,75	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS CLASSIC HEXÁGONO INTERNO Ø4,2	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS CLASSIC HEXÁGONO INTERNO Ø5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS CLASSIC HEXÁGONO INTERNO Ø6	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS PREMIUM HEXÁGONO INTERNO Ø3,3	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS PREMIUM HEXÁGONO INTERNO Ø3,75	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS PREMIUM HEXÁGONO INTERNO Ø4,2	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS PREMIUM HEXÁGONO INTERNO Ø5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ARDS PREMIUM HEXÁGONO INTERNO Ø6	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
ASTRA CONE 20°	TPDH14L52	1,4	15	5,2	2,9	3,4			1,8	2,3	Straight	45° Chamfer	1,70 Hexalobular
ASTRA CONE 45°	TPDH14L52	1,4	15	5,2	2,9	3,4			1,8	2,3	Straight	45° Chamfer	1,70 Hexalobular
ASTRA REGULAR PLATFORM	TPDH16L76	1,6	20	7,6	3,6	6,1			1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
ASTRA WIDE PLATFORM	TPDH2L90	2	25	9	4	7,5			1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
ASTRA EVOLUTION EV3.0	TPDH14L74	1,4	15	7,4	3,55	5	5,99	1,41		2,3	Conical	25° 45° Chamfer	1,70 Hexalobular
ASTRA EVOLUTION EV3.6	TPDH16L72	1,6	20	7,2	3,5	5,2	5,85	1,35		2,3	Conical	30° 45° Chamfer	1,70 Hexalobular
ASTRA EVOLUTION EV4.2	TPDH18L74	1,8	25	7,4	3,8	5,8	6,03	1,6		2,3	Conical	50° 45° Chamfer	1,70 Hexalobular
ASTRA EVOLUTION EV4.8	TPDH2L74	2	25	7,4	3,3	6			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
ASTRA EVOLUTION EV5.4	TPDH2L74	2	25	7,4	3,3	6			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
ASTRA EVOLUTION UNIABUTMENT EXTERNAL CONE	TPDH18L45	1,8	25	4,5	2,3	2,81	2,9	1,6		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
AVINENT EXTERNAL HEXAGON NARROW PLATFORM	TPDH18	1,8	25	4,4	2,75	3			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
AVINENT EXTERNAL HEXAGON REGULAR PLATFORM	TPDH2L6	2	25	6	2,7	4,5			1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
AVINENT EXTERNAL HEXAGON WIDE PLATFORM	TPDH25L67	2,5	25	6,7	4,6	5,1			1,6	2,85	Straight	Semiespherical	1,70 Hexalobular
AVINENT INTERNAL HEXAGON PLAT. 4,2	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
AVINENT INTERNAL REGULAR PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45		2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
AVINENT MULTI-UNIT	TPDH14L39	1,4	15	3,9	1,8	2,1			1,8	2,4	Straight	45° Chamfer	1,70 Hexalobular
BEGO EXTERNAL HEXAGON MINI PLATFORM 3,5	TPDH16L6	1,6	20	5,9	3	4,4			1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
BEGO INTERNAL HEXAGON PLATFORM 3,75	TPDH18L75	1,8	25	7,5	3,3	6,1			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
BEGO INTERNAL HEXAGON PLATFORM 4,1	TPDH18L75	1,8	25	7,5	3,3	6,1			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
BEGO INTERNAL HEXAGON PLATFORM 4,5	TPDH18L75	1,8	25	7,5	3,3	6,1			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
BEGO INTERNAL HEXAGON PLATFORM 5,5	TPDH18L75	1,8	25	7,5	3,3	6,1			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
BIOGENESIS 3ICON EXTERNAL HEXAGON NARROW PLATFORM	TPDH16L6	1,6	20	5,9	3	4,4			1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
BIOGENESIS 3ICON EXTERNAL HEXAGON REGULAR PLATFORM	TPDH2L6	2	25	6	2,7	4,5			1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
BIOGENESIS 3ICON EXTERNAL HEXAGON WIDE PLATFORM	TPDH2L6	2	25	6	2,7	4,5			1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
BIOGENESIS ATICON INTERNAL HEXAGON	TPDH2L90	2	25	9	4	7,5			1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
BIOGENESIS ITICON INTERNAL OCTAGON	TPDH2L82	2	25	8,2	4,7	6,7			1,5	2,4	Straight	Semiespherical	1,70 Hexalobular
BIOGENESIS SUPERICON WIDE EXTERNAL	TPDH25L67	2,5	25	6,7	4,6	5,1			1,6	2,85	Straight	Semiespherical	1,70 Hexalobular
BIOGENESIS SUPERICON WIDE INTERNAL	TPDH2L74	2	25	7,4	3,3	6			1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
BIOHORIZONS EXTERNAL HEXAGON NARROW PLATFORM	TPDH18L65	1,8	25	6,5	2,83	5			1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
BIOHORIZONS EXTERNAL HEXAGON REGULAR PLATFORM	TPDH2L65	2	25	6,5	2,7	5			1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)			HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION		
				TOTAL	THREAD	A						B	C
BIOHORIZONS EXTERNAL HEXAGON WIDE PLATFORM	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BIOHORIZONS INTERNAL HEXAGON PLATFORM 3,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BIOHORIZONS INTERNAL HEXAGON PLATFORM 4,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BIOHORIZONS INTERNAL HEXAGON PLATFORM 5,7	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BILOK EXTERNAL HEXAGON PLATFORM 3,5	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIOMET 3I CERTAIN NARROW PLATFORM	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
BIOMET 3I CERTAIN REGULAR PLATFORM	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
BIOMET 3I CERTAIN WIDE PLATFORM	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
BIOMET 3I OSSEOTITE MICRO MINI (NARROW PLATFORM)	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
BIOMET 3I OSSEOTITE REGULAR PLATFORM	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIOMET 3I OSSEOTITE WIDE PLATFORM	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIONER IKELT/BIKELT HEXÁGONO EXTERNO RP	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIONER IKELT HEXÁGONO EXTERNO WP Ø5,0	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIONER IKELT HEXÁGONO EXTERNO WP Ø5,1	TPDH25L67	2,5	25	6,7	4,6	5,1		1,6	2,85	Straight		Semiespherical	1,70 Hexalobular
BIONER HIKELT HEXÁGONO INTERNO RP Ø3,8	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BIONER HIKELT HEXÁGONO INTERNO RP Ø4,7	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BIONER TRANSEPITELIAL A-5M MULTI-UNIT RP	TPDH14L39	1,4	15	3,9	1,8	2,1		1,8	2,4	Straight		45° Chamfer	1,70 Hexalobular
BIONER TOP-DM Ø3,5	TPDH16L73	1,6	20	7,3	2,2	4,87	5,56	1,74	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
BIONER TOP-DM Ø4,0	TPDH16L73	1,6	20	7,3	2,2	4,87	5,56	1,74	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
BIONER TOP-DM Ø5,0	TPDH16L73	1,6	20	7,3	2,2	4,87	5,56	1,74	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
BIONNOVATION EXTERNAL HEXAGON TD PLATFORM 4,1	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
BIONNOVATION EXTERNAL HEXAGON TD PLATFORM 5	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
BIONNOVATION EXTERNAL HEXAGON SP CLASSIC CYLINDRICAL	TPDH18L65	1,8	25	6,5	2,83	5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIONNOVATION EXTERNAL HEXAGON RP CLASSIC CYLINDRICAL	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
BIONNOVATION EXTERNAL HEXAGON WP CLASSIC CYLINDRICAL	TPDH25L67	2,5	25	6,7	4,6	5,1		1,6	2,85	Straight		Semiespherical	1,70 Hexalobular
BIONNOVATION EXTERNAL HEXAGON SP CONIC	TPDH18L65	1,8	25	6,5	2,83	5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIONNOVATION EXTERNAL HEXAGON RP CONIC	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
BIONNOVATION EXTERNAL HEXAGON WP CONIC	TPDH25L67	2,5	25	6,7	4,6	5,1		1,6	2,85	Straight		Semiespherical	1,70 Hexalobular
BIONNOVATION INTERNAL HEXAGON SP CONIC	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
BIONNOVATION INTERNAL HEXAGON RP CONIC	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
BIONNOVATION INTERNAL HEXAGON WP CONIC	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
BIOTECH EXTERNAL OCTOGON Ø3,6	TPDH2L79	2	25	7,9	3,3	6,33	6,5	1,4	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
BIOTEC SPR/CIM Ø3,3	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BIOTEC SPR/CIM Ø3,75	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BIOTEC SPR/SPTT/CIM Ø4,2	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BIOTEC SPR/SPTT/CIM Ø5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BIOTEC SPR/SPTT/CIM Ø6	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BIOTEM AR CONICAL PLATFORM MINI Ø3,4	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIOTEM AR CONICAL PLATFORM RP	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIOTEM AR CONICAL PLATFORM WP	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIOTEM BR EXTERNAL HEXAGON NP	TPDH16L6	1,6	20	5,9	3	4,4		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIOTEM BR EXTERNAL HEXAGON RP	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
BIOTEM BR EXTERNAL HEXAGON WP	TPDH25L67	2,5	25	6,7	4,6	5,1		1,6	2,85	Straight		Semiespherical	1,70 Hexalobular
BLUE SKY BIO ONE STAGE REGULAR PLATFORM	TPDH2L74	2	25	7,4	3,3	6		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
BLUE SKY BIO INTERNAL HEX NARROW PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BLUE SKY BIO INTERNAL HEX REGULAR PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BLUE SKY BIO INTERNAL HEX WIDE PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BLUE SKY BIO TRILOBE NARROW PLATFORM	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
BLUE SKY BIO TRILOBE REGULAR PLATFORM	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)			HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION		
				TOTAL	THREAD	A						B	C
BLUE SKY BIO TRILOBE WIDE PLATFORM	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BLUE SKY BIO MAX NARROW PLATFORM	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTI MULTI-IM	TPDH14L39	1,4	15	3,9	1,8	2,1		1,8	2,4	Straight	45° Chamfer	1,70 Hexalobular	
BTI EXTERNAL NARROW PLATFORM TINY	TPDH2L51	2	25	5,1	3,1	3,6		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTI EXTERNAL REGULAR PLATFORM	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTI EXTERNAL WIDE PLATFORM	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTI INTERNAL REGULAR PLATFORM	TPDH18L65	1,8	25	6,5	2,83	5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTI INTERNAL WIDE PLATFORM	TPDH18L65	1,8	25	6,5	2,83	5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KLASSIC EXTERNAL HEXAGON PLATFORM 3,4	TPDH18L65	1,8	25	6,5	2,83	5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KLASSIC EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KLASSIC EXTERNAL HEXAGON PLATFORM 5	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KONIC EXTERNAL HEXAGON PLATFORM 3,4	TPDH18L65	1,8	25	6,5	2,83	5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KONIC EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KONIC EXTERNAL HEXAGON PLATFORM 5	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KLASSIC INTERNAL HEXAGON PLATFORM 3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KLASSIC INTERNAL HEXAGON PLATFORM 4,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KLASSIC INTERNAL HEXAGON PLATFORM 5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KONIC INTERNAL HEXAGON PLATFORM 3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KONIC INTERNAL HEXAGON PLATFORM 4,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK KONIC INTERNAL HEXAGON PLATFORM 5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
BTK SAFE TL INTERNAL OCTAGON PLATFORM 4,8	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular	
BW EXTERNAL HEXAGON NARROW PLATFORM Ø3,3	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
BW EXTERNAL HEXAGON REGULAR PLATFORM Ø4,1	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
BW EXTERNAL HEXAGON WIDE PLATFORM Ø5	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
CAMLOG SCREW LINE 3.3	TPDH16L94	1,6	20	9,4	2,9	7,65	8	1,4	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
CAMLOG SCREW LINE 3.8	TPDH16L94	1,6	20	9,4	2,9	7,65	8	1,4	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
CAMLOG SCREW LINE 4.3	TPDH16L94	1,6	20	9,4	2,9	7,65	8	1,4	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
CAMLOG SCREW LINE 5.0		#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A
CAMLOG SCREW LINE 6.0		#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A	#N/A
CAMLOG SCREW LINE 5.0	TPDH2L94	2	25	9,4	3	2,3	7,85	8	1,4	Conical	45°	45° chamfer	1,70 Hexalobular
CAMLOG SCREW LINE 6.0	TPDH2L94	2	25	9,4	3	2,3	7,85	8	1,4	Conical	45°	45° chamfer	1,70 Hexalobular
COLOSSO MORSE CONE/INTERNAL HEXAGON Ø3,5	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
COLOSSO MORSE CONE/INTERNAL HEXAGON Ø4,0	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
COLOSSO MORSE CONE/INTERNAL HEXAGON Ø5,0	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
COLOSSO EXTERNAL HEXAGON Ø3,5	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
COLOSSO EXTERNAL HEXAGON Ø4,1	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
COLOSSO EXTERNAL HEXAGON Ø5,0	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
CONEXAO AR MORSE Ø3,75	TPDH16L73	1,6	20	7,3	2,2	4,87	5,56	1,74	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
CONEXAO AR MORSE Ø4	TPDH16L73	1,6	20	7,3	2,2	4,87	5,56	1,74	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
CONEXAO MASTER ACTIVE EXTERNAL HEXAGON Ø3,3	TPDH2	2	25	4,7	3	3,3		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO MASTER ACTIVE EXTERNAL HEXAGON Ø3,75	TPDH2	2	25	4,7	3	3,3		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO MASTER ACTIVE EXTERNAL HEXAGON Ø4	TPDH2	2	25	4,7	3	3,3		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO MASTER ACTIVE EXTERNAL HEXAGON Ø5	TPDH2	2	25	4,7	3	3,3		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO MASTER POROUS EXTERNAL HEXAGON Ø3,3	TPDH2	2	25	4,7	3	3,3		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO MASTER POROUS EXTERNAL HEXAGON Ø3,75	TPDH2	2	25	4,7	3	3,3		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO MASTER POROUS EXTERNAL HEXAGON Ø4	TPDH2	2	25	4,7	3	3,3		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO MASTER POROUS EXTERNAL HEXAGON Ø5	TPDH2	2	25	4,7	3	3,3		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO GRIP ACTIVES EXTERNAL HEXAGON Ø3,75	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
CONEXAO GRIP ACTIVES EXTERNAL HEXAGON Ø4	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)					HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION
				TOTAL	THREAD	A	B	C					
CONEXAO GRIP POROUS EXTERNAL HEXAGON Ø3,75	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
CONEXAO GRIP POROUS EXTERNAL HEXAGON Ø4	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
CONEXAO CONECT AR ACTIVES INTERNAL HEXAGON Ø3,75	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO CONECT AR ACTIVES INTERNAL HEXAGON Ø4	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO CONECT AR ACTIVES INTERNAL HEXAGON Ø5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO CONECT AR INTERNAL HEXAGON Ø3,75	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO CONECT AR INTERNAL HEXAGON Ø4	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO CONECT AR INTERNAL HEXAGON Ø5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO CONICO ACTIVES INTERNAL HEXAGON Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO CONICO ACTIVES INTERNAL HEXAGON Ø4,3	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO CONICO ACTIVES INTERNAL HEXAGON Ø5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO CONICO INTERNAL HEXAGON Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO CONICO INTERNAL HEXAGON Ø4	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
CONEXAO CONICO INTERNAL HEXAGON Ø5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
DENSPLY XIVE NP Ø3,4	TPDH16L81	1,6	20	8,1	6,35	6,35	6,72	1,38	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
DENSPLY XIVE RP Ø3,8	TPDH16L81	1,6	20	8,1	6,35	6,35	6,72	1,38	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
DENSPLY XIVE WP	TPDH16L81	1,6	20	8,1	6,35	6,35	6,72	1,38	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
DENTAL IMPLANT KA INTERNAL HEXAGON PLATFORM 3,0 Ø3,3	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTAL IMPLANT KA INTERNAL HEXAGON PLATFORM 3,0 Ø3,7	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTAL IMPLANT KA INTERNAL HEXAGON PLATFORM 3,5 Ø4,1	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTAL IMPLANT KA INTERNAL HEXAGON PLATFORM 3,5 Ø4,7	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTAL IMPLANT KA INTERNAL HEXAGON PLATFORM 4,5 Ø6	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTAL TECH INTERNAL HEXAGON REGULAR PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
DENTAURUM TIOLOGIC M	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTAURUM TIOLOGIC S	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTAURUM TIOLOGIC L	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTIS S-CLEAN RP	TPDH2L79	2	25	7,9	3,3	6,33	6,5	1,4	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
DENTIUM NR LINE PLATFORM 3,5	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTIUM NR LINE PLATFORM 4,3	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTIUM SIMPLE LINE PLATFORM 4,8	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
DENTIUM SIMPLE LINE PLATFORM 6,5	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
DENTIUM SUPER LINE SE	TPDH2L79	2	25	7,9	3,3	6,33	6,5	1,4	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
DENTOFLEX SUPREME MORSE CONE INTERNAL HEXAGON Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTOFLEX SUPREME MORSE CONE INTERNAL HEXAGON Ø3,75	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTOFLEX SUPREME MORSE CONE INTERNAL HEXAGON Ø4,25	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
DENTOFLEX SUPREME MORSE CONE INTERNAL HEXAGON Ø5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
DIO FTM EXTERNAL HEXAGON NARROW PLATFORM Ø3,5	TPDH16L6	1,6	20	5,9	3	4,4		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
DIO FTM EXTERNAL HEXAGON REGULAR PLATFORM Ø4,1	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
DIO FTM EXTERNAL HEXAGON WIDE PLATFORM Ø5,1	TPDH25L67	2,5	25	6,7	4,6	5,1		1,6	2,85	Straight		Semiespherical	1,70 Hexalobular
DIO SM SYSTEM NARROW PLATFORM	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
DIO SM SYSTEM REGULAR PLATFORM	TPDH2L74	2	25	7,4	3,3	6		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
DIO UF NP	TPDH14L67	1,4	15	6,7	2,31	5	5,45	1,25	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
DIO UF RP	TPDH2L79	2	25	7,9	3,3	6,33	6,5	1,4	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
DSP BIOMEDICAL EXTERNAL HEXAGON PLATFORM 3,3	TPDH18L65	1,8	25	6,5	2,83	5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
DSP BIOMEDICAL EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
DSP BIOMEDICAL EXTERNAL HEXAGON PLATFORM 5	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
DSP BIOMEDICAL CYLINDRICAL MORSE CONE	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
DSP BIOMEDICAL CONIC MORSE CONE	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
DSP BIOMEDICAL HEX MORSE CONE PLATFORM 3,3	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)			HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION
				TOTAL	THREAD	A					
DSP BIOMEDICAL HEX MORSE CONE PLATFORM 3,75	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
DSP BIOMEDICAL HEX MORSE CONE PLATFORM 4	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
DSP BIOMEDICAL HEX MORSE CONE PLATFORM 5	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
DSP BIOMEDICAL HEX MORSE CONE CRESTAL LEVEL Ø3,5	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
DYNA UNIVERSAL 3,6	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
DYNA OCTAGON	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
ECKERMANN ALL-SPIRAL	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
ECKERMAN HEXAGON Ø 4,1	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
ECKERMANN TRIPLO	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
ELITE MEDICAL EXTERNAL HEXAGON NP	TPDH16L6	1,6	20	5,9	3	4,4	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
ELITE MEDICAL EXTERNAL HEXAGON RP	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
ELITE MEDICAL EXTERNAL HEXAGON WP	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular
EMFILS NOVO COLOSSO MORSE CONE Ø3,5	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular
EMFILS NOVO COLOSSO MORSE CONE Ø4	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular
EMFILS NOVO COLOSSO MORSE CONE Ø5	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular
EMFILS NOVO COLOSSO INTERNAL HEXAGON Ø3,5	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular
EMFILS NOVO COLOSSO INTERNAL HEXAGON Ø4	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular
EMFILS NOVO COLOSSO INTERNAL HEXAGON Ø5	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular
EMFILS NOVO COLOSSO EXTERNAL HEXAGON Ø3,5	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular
EMFILS NOVO COLOSSO EXTERNAL HEXAGON Ø4	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular
EMFILS NOVO COLOSSO EXTERNAL HEXAGON Ø5	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular
EUROTEKNIKA AESTHETICA INTERNAL OCTAGON Ø3,6	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
EUROTEKNIKA AESTHETICA INTERNAL OCTAGON Ø4,1	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
EUROTEKNIKA AESTHETICA INTERNAL OCTAGON Ø4,8	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
EUROTEKNIKA UNIVERSAL RP	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
EUROTEKNIKA NATURACTIS NARROW PLATFORM	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
EUROTEKNIKA NATURACTIS REGULAR PLATFORM	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
EUROTEKNIKA NATURACTIS WIDE PLATFORM	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
EUROTEKNIKA NATEA/NATURALL NP	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
EUROTEKNIKA NATEA/NATURALL RP	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
EUROTEKNIKA NATEA/NATURALL WP	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
EUROTEKNIKA UNEVA 4,1	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
EUROTEKNIKA UNEVA NE RP	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
FRIALIT MULTI-UNIT MP	TPDH16	1,6	20	4,4	2,5	2,9	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
GALIMPLANT EXTERNAL HEXAGON PLATFORM Ø4	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular
GALIMPLANT INTERNAL HEXAGON PLATFORM Ø3,5	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular
GALIMPLANT INTERNAL HEXAGON PLATFORM Ø4,0	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
GMI FRONTIER STANDARD	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
GMI FRONTIER WIDE	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
GMI PHOENIX NARROW 3,3	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
GMI PHOENIX STANDARD 4,1	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
GMI PHOENIX WIDE 5,1 (E/C 2,7)	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
GMI PHOENIX WIDE 5,1 (E/C 3,4)	TPDH25	2,5	25	5,4	3,8	4,1	1,3	2,85	Straight	Semiespherical	1,70 Hexalobular
GMI PHOENIX SUPER WIDE 6	TPDH25	2,5	25	5,4	3,8	4,1	1,3	2,85	Straight	Semiespherical	1,70 Hexalobular
GT MEDICAL BEST-FIT EXTERNAL HEXAGON NARROW PLATFORM	TPDH16L6	1,6	20	5,9	3	4,4	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
GT MEDICAL BEST-FIT EXTERNAL HEXAGON REGULAR PLATFORM	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
GT MEDICAL BEST-FIT EXTERNAL HEXAGON WIDE PLATFORM	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular
GT MEDICAL BEST-FIT INTERNAL HEXAGON Ø3,7	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
GT MEDICAL BEST-FIT INTERNAL HEXAGON Ø4,3	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)					HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION
				TOTAL	THREAD	A	B	C					
GT MEDICAL BEST-FIT INTERNAL HEXAGON Ø4,8	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
GT MEDICAL BEST-FIT INTERNAL OCTOGON	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
HI-TEC SELF THREAD SLT NP Ø3,75	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
HI-TEC SELF THREAD GTI RP Ø4,2	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
HI-TEC LGI PLUS NP Ø4,3	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
IMPLACIL DE BORTOLI EXTERNAL HEXAGON PLATFORM Ø3,5	TPDH16L6	1,6	20	5,9	3	4,4		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLACIL DE BORTOLI EXTERNAL HEXAGON PLATFORM Ø4	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
IMPLACIL DE BORTOLI EXTERNAL HEXAGON PLATFORM Ø5	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
IMPLACIL DE BORTOLI INTERNAL HEXAGON PLATFORM Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLACIL DE BORTOLI INTERNAL HEXAGON PLATFORM Ø4	TPDH2L74	2	25	7,4	3,3	6		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLACIL DE BORTOLI INTERNAL HEXAGON PLATFORM Ø5	TPDH2L74	2	25	7,4	3,3	6		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLACIL DE BORTOLI MORSE CONE Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLACIL DE BORTOLI MORSE CONE Ø4	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLACIL DE BORTOLI MORSE CONE Ø4,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLACIL DE BORTOLI MORSE CONE Ø5,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT REACTIVE Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT REACTIVE Ø4,3	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT REACTIVE Ø5	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT REPLANT Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT REPLANT Ø4,3	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT REPLANT Ø5	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT REPLANT Ø6	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT REPLUS Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT REPLUS Ø4,3	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT REPLUS Ø5	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT LEGACY Ø3,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
IMPLANT DIRECT LEGACY Ø4,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
IMPLANT DIRECT LEGACY Ø5,7	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
IMPLANT DIRECT SWISSPLANT REGULAR PLATFORM	TPDH2L74	2	25	7,4	3,3	6		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular
IMPLANT DIRECT INTERACTIVE NARROW PLATFORM	TPDH16L73	1,6	20	7,3	2,2	4,87	5,56	1,74	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
IMPLANT DIRECT INTERACTIVE REGULAR PLATFORM	TPDH2L75	2	25	7,5	2,75	5,93	6,18	1,32	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
IMPLASA HÖSCHT DOUBLEGOTIC Ø3,3	TPDH18L74	1,8	25	7,4	3,8	5,8	6,03	1,6	2,3	Conical	50°	45° Chamfer	1,70 Hexalobular
IMPLASA HÖSCHT DOUBLEGOTIC/GOTIC/NORDIC Ø3,75	TPDH18L74	1,8	25	7,4	3,8	5,8	6,03	1,6	2,3	Conical	50°	45° Chamfer	1,70 Hexalobular
IMPLASA HÖSCHT DOUBLEGOTIC/GOTIC/NORDIC Ø4,2	TPDH18L74	1,8	25	7,4	3,8	5,8	6,03	1,6	2,3	Conical	50°	45° Chamfer	1,70 Hexalobular
IMPLASA HÖSCHT DOUBLEGOTIC/GOTIC/NORDIC Ø5	TPDH18L74	1,8	25	7,4	3,8	5,8	6,03	1,6	2,3	Conical	50°	45° Chamfer	1,70 Hexalobular
IMPLASA HÖSCHT DOUBLEGOTIC/GOTIC/NORDIC Ø6	TPDH18L74	1,8	25	7,4	3,8	5,8	6,03	1,6	2,3	Conical	50°	45° Chamfer	1,70 Hexalobular
INNOVA ENDOPORE EXTERNAL HEXAGON REGULAR PLATFORM	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
JDENTALCARE JDEVOLUTION INTERNAL HEXAGON Ø3,7	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
JDENTALCARE JDEVOLUTION INTERNAL HEXAGON Ø4,3	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
JDENTALCARE JDEVOLUTION INTERNAL HEXAGON Ø5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
JDENTALCARE JDEVOLUTION INTERNAL HEXAGON Ø6	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
JDENTALCARE JDEVOLUTIONS Ø3,2	TPDH16L73	1,6	20	7,3	2,2	4,87	5,56	1,74	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
KENTRON INTERNAL HEXAGON NARROW PLATFORM 3,8	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
KENTRON INTERNAL HEXAGON REGULAR PLATFORM 4,5	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight		Semiespherical	1,70 Hexalobular
KEYSTONE PRIMACONNEX SD	TPDH18L65	1,8	25	6,5	2,83	5			2,3	Straight		Semiespherical	1,70 Hexalobular
KEYSTONE PRIMACONNEX RD	TPDH18L65	1,8	25	6,5	2,83	5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
KEYSTONE PRIMACONNEX WD	TPDH18L65	1,8	25	6,5	2,83	5		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
KEYSTONE PRIMACONNEX MULTI UNIT SD/RD	TPDH16	1,6	20	4,4	2,5	2,9		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
KEYSTONE PRIMACONNEX MULTI UNIT WD	TPDH16	1,6	20	4,4	2,5	2,9		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)					HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION
				TOTAL	THREAD	A	B	C					
KEYSTONE RESTORE NP	TPDH18	1,8	25	4,4	2,75	3		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
KEYSTONE RESTORE RP	TPDH18	1,8	25	4,4	2,75	3		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
KEYSTONE RESTORE WP	TPDH18	1,8	25	4,4	2,75	3		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
KLOCKNER ESSENTIAL CONE REGULAR PLATFORM	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
KLOCKNER ESSENTIAL CONE WIDE PLATFORM	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
KLOCKNER ESSENTIAL CONE OCTACONE REGULAR PLATFORM	TPDH18	1,8	25	4,4	2,75	3		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
KLOCKNER ESSENTIAL CONE OCTACONE WIDE PLATFORM	TPDH18	1,8	25	4,4	2,75	3		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
KLOCKNER NK2-SK2	TPDH18	1,8	25	4,4	2,75	3		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
KLOCKNER KL NP	TPDH16L6	1,6	20	5,9	3	4,4		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
KLOCKNER KL RP	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular	
KLOCKNER KL WP	TPDH25L67	2,5	25	6,7	4,6	5,1		1,6	2,85	Straight	Semiespherical	1,70 Hexalobular	
LASAQ BIONIC 3.5/4.0/5.0	TPDH16L81	1,6	20	8,1	6,35	6,35	6,72	1,38	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
LEADER ITALY TIXOS INTERNAL HEXAGON PLATFORM 3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
LEADER ITALY TIXOS INTERNAL HEXAGON PLATFORM 4	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
LEADER ITALY TIXOS INTERNAL HEXAGON PLATFORM 5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
LEADER ITALY TIXOS INTERNAL HEXAGON PLATFORM 6	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
LEADER ITALY TIXOS EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular	
LEADER ITALY TIXOS EXTERNAL HEXAGON PLATFORM 5	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular	
LEADER ITALY TIXOS MORSE CONE PLATFORM 3,6	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight	Semiespherical	1,70 Hexalobular	
LIFECORE EXTERNAL HEXAGON NARROW PLATFORM	TPDH18L65	1,8	25	6,5	2,83	5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
LIFECORE EXTERNAL HEXAGON REGULAR PLATFORM	TPDH2L65	2	25	6,5	2,7	5		1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular	
LIFECORE EXTERNAL HEXAGON WIDE PLATFORM	TPDH25L67	2,5	25	6,7	4,6	5,1		1,6	2,85	Straight	Semiespherical	1,70 Hexalobular	
MEDENTIS ICX TEMPLANT Ø3,3	TPDH16L81	1,6	20	8,1	6,35	6,35	6,72	1,38	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
MEDENTIS ICX TEMPLANT Ø3,75	TPDH16L81	1,6	20	8,1	6,35	6,35	6,72	1,38	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
MEDENTIS ICX TEMPLANT Ø4,1	TPDH16L81	1,6	20	8,1	6,35	6,35	6,72	1,38	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
MEDENTIS ICX TEMPLANT Ø4,8	TPDH16L81	1,6	20	8,1	6,35	6,35	6,72	1,38	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
MEGAGEN ANYONE INTERNAL RP	TPDH2L79	2	25	7,9	3,3	6,33	6,5	1,4	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
MEGAGEN ANYRIDGE INTERNAL HEXAGON SMALL DIAMETER 3,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN ANYRIDGE INTERNAL HEXAGON REGULAR DIAMETER 4	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN ANYRIDGE INTERNAL HEXAGON REGULAR DIAMETER 4,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN ANYRIDGE INTERNAL HEXAGON WIDE DIAMETER 5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN ANYRIDGE INTERNAL HEXAGON WIDE DIAMETER 5,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN ANYRIDGE INTERNAL HEXAGON SUPERWIDE DIAMETER 6	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN ANYRIDGE INTERNAL HEXAGON SUPERWIDE DIAMETER 6,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN ANYRIDGE INTERNAL HEXAGON SUPERWIDE DIAMETER 7	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN ANYRIDGE INTERNAL HEXAGON SUPERWIDE DIAMETER 7,5	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN ANYRIDGE INTERNAL HEXAGON SUPERWIDE DIAMETER 8	TPDH18L75	1,8	25	7,5	3,3	6,1		1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN EZ PLUS INTERNAL CONECTION SMALL Ø3,5	TPDH16L84	1,6	20	8,4	3,5	6,8		1,6	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN EZ PLUS INTERNAL CONECTION REGULAR Ø4,1	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN EZ PLUS INTERNAL CONECTION WIDE Ø5	TPDH2L90	2	25	9	4	7,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN RESCUE INTERNAL PLATFORM Ø4,5	TPDH25	2,5	25	5,4	3,8	4,1		1,3	2,85	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN RESCUE INTERNAL PLATFORM Ø5	TPDH25	2,5	25	5,4	3,8	4,1		1,3	2,85	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN RESCUE INTERNAL PLATFORM Ø5,5	TPDH25	2,5	25	5,4	3,8	4,1		1,3	2,85	Straight	Semiespherical	1,70 Hexalobular	
MEGAGEN RESCUE INTERNAL PLATFORM 6	TPDH25	2,5	25	5,4	3,8	4,1		1,3	2,85	Straight	Semiespherical	1,70 Hexalobular	
MICRODENT MINI EXTERNAL HEXAGON Ø3,5	TPDH18L65	1,8	25	6,5	2,83	5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
MICRODENT SYSTEM EXTERNAL HEXAGON Ø4,2	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
MICRODENT SYSTEM EXTERNAL HEXAGON Ø5,1	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
MICRODENT UNIVERSAL EXTERNAL HEXAGON Ø4,1	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
MICRODENT UNIVERSAL EXTERNAL HEXAGON Ø5,1	TPDH2L6	2	25	6	2,7	4,5		1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)			HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION		
				TOTAL	THREAD	A						B	C
MICRODENT MK PLATFORM WIDE	TPDH25	2,5	25	5,4	3,8	4,1	1,3	2,85	Straight	Semiespherical	1,70 Hexalobular		
MICRODENT GENIUS INTERNAL Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
MICRODENT GENIUS INTERNAL Ø4	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
MICRODENT GENIUS INTERNAL Ø4,5	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
MICRODENT TRYLOGIC INTERNAL Ø3,8	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
MICRODENT TRYLOGIC INTERNAL Ø4,3	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
MICRODENT TRYLOGIC INTERNAL Ø5,1	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
MICRODENT EKTOS Ø3,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
MIS EXTERNAL HEXAGON 3,75	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
MIS EXTERNAL HEXAGON 4,2	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
MIS EXTERNAL HEXAGON 5	TPDH2	2	25	4,7	3	3,3	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
MIS C1 NARROW PLATFORM	TPDH16L71	1,6	20	7,1	2,8	5,2	5,53	1,57	2,3	Conical	60°	45° Chamfer	1,7 Hexalobular
MIS C1 STANDARD PLATFORM	TPDH17L73	N1-72	25	7,3	2,5	5,5	5,77	1,53	2,3	Conical	60°	45° Chamfer	1,7 Hexalobular
MIS C1 WIDE PLATFORM	TPDH17L73	N1-72	25	7,3	2,5	5,5	5,77	1,53	2,3	Conical	60°	45° Chamfer	1,7 Hexalobular
MIS V3 STANDARD PLATFORM	TPDH17L73	N1-72	25	7,3	2,5	5,5	5,77	1,53	2,3	Conical	60°	45° Chamfer	1,7 Hexalobular
MIS SEVEN/M4 NARROW PLATFORM	TPDH16L71	1,6	20	7,1	2,8	5,2	5,53	1,57	2,3	Conical	60°	45° Chamfer	1,7 Hexalobular
MIS SEVEN/M4 STANDARD PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
MIS SEVEN/M4 WIDE PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
MIS MULTI-UNIT	TPDH16	1,6	20	4,4	2,5	2,9	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
MOZO GRAU EXTERNAL HEXAGON NARROW PLATFORM	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
MOZO GRAU EXTERNAL HEXAGON REGULAR PLATFORM	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
MOZO GRAU EXTERNAL HEXAGON WIDE PLATFORM	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular		
MOZO GRAU INHEX MINI	TPDH14L67	1,4	15	6,7	2,31	5	5,45	1,25	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
MOZO GRAU INHEX STANDARD	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
MOZO GRAU INHEX MAXI	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
MPI HE EXTERNAL HEXAGON NARROW PLATFORM Ø3,3/3,5	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
MPI HE EXTERNAL HEXAGON REGULAR PLATFORM Ø4,1	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
MPI HE EXTERNAL HEXAGON WIDE PLATFORM	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
MPI CM INTERNAL HEXAGON REGULAR PLATFORM	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
MPI CM INTERNAL HEXAGON WIDE PLATFORM	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
NEOBIOTECH EXTERNAL HEXAGON EB NP	TPDH16L6	1,6	20	5,9	3	4,4	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
NEOBIOTECH EXTERNAL HEXAGON EB RP	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
NEOBIOTECH EXTERNAL HEXAGON EB WP (3I)	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
NEOBIOTECH EXTERNAL HEXAGON EB WP (BRANEMARK)	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular		
NEOBIOTECH EXTERNAL HEXAGON EB S-WIDE (BRANEMARK)	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular		
NEOBIOTECH MULTI-UNIT REGULAR PLATFORM	TPDH14L39	1,4	15	3,9	1,8	2,1	1,8	2,4	Straight	45° Chamfer	1,70 Hexalobular		
NEOBIOTECH IS CONE INTERNAL HEXAGON REGULAR PLATFORM	TPDH2L79	2	25	7,9	3,3	6,33	6,5	1,4	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
NEOBIOTECH IS CONE INTERNAL HEXAGON WIDE PLATFORM	TPDH2L79	2	25	7,9	3,3	6,33	6,5	1,4	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
NEODENT EXTERNAL HEXAGON Ø3,3	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
NEODENT EXTERNAL HEXAGON Ø4,1	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
NEODENT EXTERNAL HEXAGON Ø4,3	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
NEODENT EXTERNAL HEXAGON Ø5	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular		
NEODENT INTERNAL HEXAGON Ø4,3	TPDH16L81	1,6	20	8,1	6,35	6,35	6,72	1,38	2,3	Conical	45°	45° Chamfer	1,70 Hexalobular
NEODENT CONE MINI-PILAR	TPDH14L39	1,4	15	3,9	1,8	2,1	1,8	2,4	Straight	45° Chamfer	1,70 Hexalobular		
NEOSS PROACTIVE (3,5/4/4,5/5/5,5)	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
NOBEL BIOCARE ACTIVE NARROW PLATFORM	TPDH16L73	1,6	20	7,3	2,2	4,87	5,56	1,74	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
NOBEL BIOCARE ACTIVE REGULAR PLATFORM	TPDH2L75	2	25	7,5	2,75	5,93	6,18	1,32	2,3	Conical	35°	45° Chamfer	1,70 Hexalobular
NOBEL BIOCARE BRANEMARK NARROW PLATFORM	TPDH16L6	1,6	20	5,9	3	4,4	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
NOBEL BIOCARE BRANEMARK REGULAR PLATFORM	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)			HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION	
				TOTAL	THREAD	A						B
NOBEL BIOCARE BRANEMARK WIDE PLATFORM	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular	
NOBEL BIOCARE MULTI-UNIT REGULAR PLATFORM	TPDH14L39	1,4	15	3,9	1,8	2,1	1,8	2,4	Straight	45° Chamfer	1,70 Hexalobular	
NOBEL BIOCARE MULTI-UNIT WIDE PLATFORM	TPDH18	1,8	25	4,4	2,75	3	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NOBEL BIOCARE REPLACE NARROW PLATFORM	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NOBEL BIOCARE REPLACE REGULAR PLATFORM	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
NOBEL BIOCARE REPLACE WIDE PLATFORM	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
NOBEL BIOCARE STANDARD	TPDH14L52	1,4	15	5,2	2,9	3,4	1,8	2,3	Straight	45° Chamfer	1,70 Hexalobular	
NUCLEOSS T5 SHORT INTERNAL HEXAGON Ø4,8	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NUCLEOSS T5 SHORT INTERNAL HEXAGON Ø5,5	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NUCLEOSS T5 SHORT INTERNAL HEXAGON Ø6,2	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NUCLEOSS T4 BONE LEVEL INTERNAL HEXAGON Ø3,4	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NUCLEOSS T4 BONE LEVEL INTERNAL HEXAGON Ø3,8	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NUCLEOSS T4 BONE LEVEL INTERNAL HEXAGON Ø4,2	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NUCLEOSS T4 BONE LEVEL INTERNAL HEXAGON Ø5	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NUCLEOSS T3 TISSUE LEVEL INTERNAL HEXAGON Ø3,4	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NUCLEOSS T3 TISSUE LEVEL INTERNAL HEXAGON Ø3,8	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NUCLEOSS T3 TISSUE LEVEL INTERNAL HEXAGON Ø4,2	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
NUCLEOSS T3 TISSUE LEVEL INTERNAL HEXAGON Ø5	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSSEOLIFE OEX P2 SERIE Ø3,3	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSSEOLIFE OEX P3 SERIE Ø3,75	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular	
OSSEOLIFE OEX P4 SERIE Ø4,25	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular	
OSSEOLIFE OEX P5 SERIE Ø5	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular	
OSSEOLIFE OIN SERIE P3 Ø3,75	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSSEOLIFE OIN SERIE P4 Ø4,25	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSSEOLIFE OIN SERIE P5 Ø5	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSSTEM TS NP	TPDH16L94	1,6	20	9,4	2,9	7,65	8	1,4	2,3	Conical	45° 45° Chamfer	1,70 Hexalobular
OSSTEM TS RP	TPDH2L79	2	25	7,9	3,3	6,33	6,5	1,4	2,3	Conical	45° 45° Chamfer	1,70 Hexalobular
OSSTEM SS Ø4,8	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSSTEM SS Ø 6	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSSTEM US MINI	TPDH16L6	1,6	20	5,9	3	4,4	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSSTEM US REGULAR	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSSTEM US WIDE	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular	
OSTEOPLUS ARGENTIUM	TPDH18	1,8	25	4,4	2,75	3	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSTEOPLUS IRIIDIUM 3.5	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSTEOPLUS IRIIDIUM 4.2	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular	
OSTEOPLUS SHI INTERNAL HEX NP Ø3,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70° 45° Chamfer	1,70 Hexalobular
OSTEOPLUS SHE EXTERNAL HEX NP	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular	
OSTEOPLUS SHE EXTERNAL HEX RP	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OT MEDICAL OT-F ² INTERNAL Ø3,4	TPDH16L6	1,6	20	5,9	3	4,4	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OT MEDICAL OT-F ² INTERNAL Ø3,8	TPDH16L6	1,6	20	5,9	3	4,4	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OT MEDICAL OT-F ² INTERNAL Ø4,1	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OT MEDICAL OT-F ² INTERNAL Ø5,0	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OT MEDICAL OT-F ³ INTERNAL Ø3,4	TPDH16L6	1,6	20	5,9	3	4,4	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OT MEDICAL OT-F ³ INTERNAL Ø3,8	TPDH16L6	1,6	20	5,9	3	4,4	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OT MEDICAL OT-F ³ INTERNAL Ø4,1	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OT MEDICAL OT-F ³ INTERNAL Ø5,0	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular	
OXY INT CLASSIC Ø3,3 MINI	TPDH2L82	2	25	8,2	4,7	6,7	1,5	2,4	Straight	Semiespherical	1,70 Hexalobular	
OXY INT PS Ø3,5 MINI	TPDH2L82	2	25	8,2	4,7	6,7	1,5	2,4	Straight	Semiespherical	1,70 Hexalobular	
OXY INT PS Ø4,0 MINI	TPDH2L82	2	25	8,2	4,7	6,7	1,5	2,4	Straight	Semiespherical	1,70 Hexalobular	





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)			HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION		
				TOTAL	THREAD	A						B	C
PHIBO TSA S3	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
PHIBO TSA S4	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
PHIBO TSH S2	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
PHIBO TSH S3	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
PHIBO TSH S4	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
PHIBO TSH S5	TPDH2	2	25	4,7	3	3,3	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
PROCLINIC EXTERNAL HEXAGON MINI NP Ø3,3	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
PROCLINIC EXTERNAL HEXAGON REGULAR PLATFORM Ø4	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
PROCLINIC EXTERNAL HEXAGON WIDE PLATFORM Ø5	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
PROCLINIC INTERNAL HEXAGON NARROW PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
PROCLINIC INTERNAL HEXAGON REGULAR PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
PROCLINIC INTERNAL OCTOGON D4.8 SP OCTA RP	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
PROCLINIC INTERNAL OCTOGON D6.5 SP OCTA WP	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
RADHEX PHE EXTERNAL HEXAGON PLATFORM Ø3,5	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
RADHEX PHE EXTERNAL HEXAGON PLATFORM Ø4,1	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
RADHEX PHI INTERNAL HEXAGON PLATFORM Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
RADHEX PHI INTERNAL HEXAGON PLATFORM Ø4,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ROSTER DENT HEX EXTERNAL HEXAGON RP	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
ROSTER DENT PLUS INTERNAL OCTOGON Ø4,8	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
ROSTER DENT CONICAL GOLD Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
ROSTER DENT CONICAL GOLD Ø4	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
ROSTER DENT CONICAL GOLD Ø5	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
ROSTER DENT CONICAL GOLD P4 Ø3,5	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
ROSTER DENT CONICAL GOLD P4 Ø4	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
ROSTER DENT CONICAL GOLD P4 Ø5	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SGS P1/P7 INTERNAL HEX Ø3,2	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
SGS P1/P7 INTERNAL HEX Ø3,75	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
SGS P1/P7 INTERNAL HEX Ø4,2	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
SGS P7 INTERNAL HEX Ø4,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
SGS P1/P7 INTERNAL HEX Ø5,0	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
SGS P1/P7 INTERNAL HEX Ø6,0	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
SGS P7N EXTERNAL CONE Ø3,0	TPDH14L74	1,4	15	7,4	3,55	5	5,99	1,41	2,3	Conical	25°	45° Chamfer	1,70 Hexalobular
SGS S6 EXTERNAL CONE Ø3,75	TPDH14L39	1,4	15	3,9	1,8	2,1	1,8	2,4	Straight	45° Chamfer	1,70 Hexalobular		
SGS S7 EXTERNAL CONE Ø3,75	TPDH14L39	1,4	15	3,9	1,8	2,1	1,8	2,4	Straight	45° Chamfer	1,70 Hexalobular		
SIGNO VINCES INTERNAL CONE DUOCON Ø3,8	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIGNO VINCES INTERNAL CONE DUOCON Ø4,6	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIGNO VINCES INTERNAL CONE COMPACT Ø4,5	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIGNO VINCES EXTERNAL HEXAGON DUO NARROW PLATFORM Ø3,8	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIGNO VINCES EXTERNAL HEXAGON DUO REGULAR PLATFORM Ø4,1	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIGNO VINCES EXTERNAL HEXAGON INTEGRA NARROW PLATFORM Ø3,8	TPDH18L65	1,8	25	6,5	2,83	5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIGNO VINCES EXTERNAL HEXAGON INTEGRA REGULAR PLATFORM Ø4,1	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIGNO VINCES EXTERNAL HEXAGON WIDE PLATFORM	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular		
SIN IMPLANTS STRONG SW MORSE PLATFORM 3,5	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIN IMPLANTS STRONG SW MORSE PLATFORM 3,8	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIN IMPLANTS STRONG SW MORSE PLATFORM 4,5	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIN IMPLANTS STRONG SW MORSE PLATFORM 5,5	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIN IMPLANTS STRONG SW HI PLATFORM 3,8	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIN IMPLANTS STRONG SW HI PLATFORM 4,5	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SIN IMPLANTS STRONG SW HI PLATFORM 5,5	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		





DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)			HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION		
				TOTAL	THREAD	A						B	C
SIN IMPLANTS TRYON EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
SIN IMPLANTS TRYON EXTERNAL HEXAGON PLATFORM 5	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
SIN IMPLANTS REVOLUTION EXTERNAL HEXAGON PLATFORM 3,4	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
SIN IMPLANTS REVOLUTION EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
SIN IMPLANTS REVOLUTION EXTERNAL HEXAGON PLATFORM 5	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
SIN IMPLANTS STYLUS EXTERNAL HEXAGON PLATFORM 3,4	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
SIN IMPLANTS STYLUS EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
SIN IMPLANTS STYLUS EXTERNAL HEXAGON PLATFORM 5	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
SIN IMPLANTS REVOLUTION COMPACT EXTERNAL HEXAGON PLATFORM 5	TPDH2	2	25	4,7	3	3,3	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SOUTHERN IMPLANT EXTERNAL HEXAGON NARROW PLATFORM BLUE Ø3.25	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
SOUTHERN IMPLANT EXTERNAL HEXAGON REGULAR PLATFORM GRAY Ø4	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SOUTHERN IMPLANT EXTERNAL HEXAGON WIDE PLATFORM YELLOW Ø5	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SOUTHERN IMPLANT EXTERNAL HEXAGON GREEN Ø6	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SOUTHERN IMPLANT EXTERNAL HEXAGON MSC	TPDH16L6	1,6	20	5,9	3	4,4	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SOUTHERN IMPLANTS INTERNAL OCTAGON Ø4,8	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
SOUTHERN IMPLANTS INTERNAL OCTAGON Ø6	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
SOUTHERN IMPLANTS TRI-NEX NP	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SOUTHERN IMPLANTS TRI-NEX RP	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SOUTHERN IMPLANTS TRI-NEX WP	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SOUTHERN IMPLANTS TRI-NEX Ø6	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SOUTHERN IMPLANTS INTERNAL HEX	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
SOUTHERN IMPLANTS COMPACT CONICAL MULTI-UNIT RP	TPDH14L39	1,4	15	3,9	1,8	2,1	1,8	2,4	Straight	45° Chamfer	1,70 Hexalobular		
SOUTHERN IMPLANTS COMPACT CONICAL MULTI-UNIT WP	TPDH14L39	1,4	15	3,9	1,8	2,1	1,8	2,4	Straight	45° Chamfer	1,70 Hexalobular		
STERNGOLD STERN EX EXTERNAL HEXAGON Ø3,75	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
STERNGOLD STERN EX EXTERNAL HEXAGON Ø4	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
STERNGOLD STERN EX EXTERNAL HEXAGON Ø5	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular		
STERNGOLD STERN IC INTERNAL OCTAGON Ø3,3	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
STERNGOLD STERN IC INTERNAL OCTAGON Ø4,1	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
STRAUMANN BONE LEVEL NP	TPDH16L78	1,6	20	7,8	2	5,36	7,03	0,81	2,3	Conical	15°	45° Chamfer	1,70 Hexalobular
STRAUMANN BONE LEVEL RP	TPDH16L78	1,6	20	7,8	2	5,36	7,03	0,81	2,3	Conical	15°	45° Chamfer	1,70 Hexalobular
STRAUMANN BONE LEVEL MULTI-BASE	TPDHST	1,4	15	4,5	2,3	2,5	3,28	1,22	2,3	Conical	30°	45° Chamfer	1,70 Hexalobular
STRAUMANN INTERNAL OCTAGON RP	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
STRAUMANN INTERNAL OCTAGON WIDE NECK WP	TPDH2L67	2	25	6,7	2,3	3,65	5,68	1,02	2,59	Conical	15°	45° Chamfer	1,70 Hexalobular
STRAUMANN SYNOCTA	TPDH2L51	2	25	5,1	3,1	3,6	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA OUTLINK ² Slim 3,0 mm	TPDH18L48	1,8	25	4,8	2,8	3,22	3,65	1,15	2,3	Conical	30°	Semiespherical	1,70 Hexalobular
SWEDEN&MARTINA OUTLINK ² 3,3 mm	TPDH18L48	1,8	25	4,8	2,8	3,22	3,65	1,15	2,3	Conical	30°	Semiespherical	1,70 Hexalobular
SWEDEN&MARTINA OUTLINK ² 3,75 mm	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA OUTLINK ² 4,1 mm	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA OUTLINK ² 5 mm	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA PREMIUM 3,3	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA PREMIUM Y KOHNO 3,8	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA PREMIUM Y KOHNO 4,25	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA PREMIUM Y KOHNO 5	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA PREMIUM Y KOHNO 6	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA STARK	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA GLOBAL Ø3,8	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA GLOBAL Ø4,3	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA GLOBAL Ø4,8	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		
SWEDEN&MARTINA GLOBAL Ø5,5	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular		



DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)			HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION
				TOTAL	THREAD	A					
SWEDEN&MARTINA SHELTA Ø3,8	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
SWEDEN&MARTINA SHELTA Ø4,25	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
SWEDEN&MARTINA SHELTA Ø5	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX MORSE CONE INTERNAL HEXAGON Ø3,5	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX MORSE CONE INTERNAL HEXAGON Ø4,3	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX MORSE CONE INTERNAL HEXAGON Ø5	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX CLASSIC-IN CR INTERNAL HEXAGON Ø3,5	TPDH16L84	1,6	20	8,4	3,5	6,8	1,6	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX CLASSIC-IN CR INTERNAL HEXAGON Ø3,75	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX CLASSIC-IN CR INTERNAL HEXAGON Ø4	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX CLASSIC-IN CR INTERNAL HEXAGON Ø5	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX ESTETIC-IN CR INTERNAL HEXAGON PLATFORM 4,3	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX ESTETIC-IN CR INTERNAL HEXAGON PLATFORM 5	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX CLASSIC-CI EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular
SYSTHEX CLASSIC-CI EXTERNAL HEXAGON PLATFORM 5	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX ESTETIC-CI EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular
SYSTHEX ESTETIC-CI EXTERNAL HEXAGON PLATFORM 5	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX FIT-HEX EXTERNAL HEXAGON PLATFORM 3,3	TPDH18	1,8	25	4,4	2,75	3	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX CLASSIC TRADICIONAL EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular
SYSTHEX CLASSIC TRADICIONAL EXTERNAL HEXAGON PLATFORM 5	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX CLASSIC CURTO EXTERNAL HEXAGON PLATFORM 3,3	TPDH18	1,8	25	4,4	2,75	3	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
SYSTHEX CLASSIC CURTO EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular
TAV DENTAL SIRIUS INTERNAL HEXAGON NARROW PLATFORM Ø3,3	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
TAV DENTAL SIRIUS INTERNAL HEXAGON STANDARD PLATFORM Ø3,75	TPDH2L82	2	25	8,2	4,7	6,7	1,5	2,4	Straight	Semiespherical	1,70 Hexalobular
TAV DENTAL SIRIUS INTERNAL HEXAGON STANDARD PLATFORM Ø4,2	TPDH2L82	2	25	8,2	4,7	6,7	1,5	2,4	Straight	Semiespherical	1,70 Hexalobular
TAV DENTAL SIRIUS INTERNAL HEXAGON STANDARD PLATFORM Ø5	TPDH2L82	2	25	8,2	4,7	6,7	1,5	2,4	Straight	Semiespherical	1,70 Hexalobular
TAV DENTAL SIRIUS INTERNAL HEXAGON STANDARD PLATFORM Ø6	TPDH2L82	2	25	8,2	4,7	6,7	1,5	2,4	Straight	Semiespherical	1,70 Hexalobular
TAV DENTAL SILHOUETTE INTERNAL HEXAGON NARROW PLATFORM Ø3,3	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
TAV DENTAL SILHOUETTE INTERNAL HEXAGON STANDARD PLATFORM Ø3,75	TPDH2L82	2	25	8,2	4,7	6,7	1,5	2,4	Straight	Semiespherical	1,70 Hexalobular
TAV DENTAL SILHOUETTE INTERNAL HEXAGON STANDARD PLATFORM Ø4,2	TPDH2L82	2	25	8,2	4,7	6,7	1,5	2,4	Straight	Semiespherical	1,70 Hexalobular
TAV DENTAL SILHOUETTE INTERNAL HEXAGON STANDARD PLATFORM Ø5	TPDH2L82	2	25	8,2	4,7	6,7	1,5	2,4	Straight	Semiespherical	1,70 Hexalobular
TAV DENTAL SILHOUETTE INTERNAL HEXAGON STANDARD PLATFORM Ø6	TPDH2L82	2	25	8,2	4,7	6,7	1,5	2,4	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX b-FIX BLACK FIX PLATFORM 3,5	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX b-FIX BLACK FIX PLATFORM 4	TPDH16L76	1,6	20	7,6	3,6	6,1	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX b-FIX BLACK FIX PLATFORM 4,5	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX b-FIX BLACK FIX PLATFORM 5	TPDH2L90	2	25	9	4	7,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX i-FIX INTERNAL HEXAGON GROOVE PLATFORM 4,1	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX i-FIX INTERNAL HEXAGON PLATFORM 4,1	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX i-FIX INTERNAL HEXAGON PLATFORM 5,1	TPDH18L75	1,8	25	7,5	3,3	6,1	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX e-FIX EXTERNAL HEXAGON INTERNAL TORQUE PLATFORM 4,1	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular
TITANIUM FIX e-FIX EXTERNAL HEXAGON INTERNAL TORQUE SILVER PLAT. 4,1	TPDH2L65	2	25	6,5	2,7	5	1,5	2,3	Straight	45° Chamfer	1,70 Hexalobular
TITANIUM FIX e-FIX EXTERNAL HEXAGON PLATFORM 4,1	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX e-FIX EXTERNAL HEXAGON PLATFORM 5,1	TPDH25L67	2,5	25	6,7	4,6	5,1	1,6	2,85	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX e-FIX SILVER PLATFORM 4,1	TPDH2	2	25	4,7	3	3,3	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX e-FIX SILVER PLATFORM 5,1	TPDH25	2,5	25	5,4	3,8	4,1	1,3	2,85	Straight	Semiespherical	1,70 Hexalobular
TITANIUM FIX c-FIX MORSE CONE PLATFORM 4,8	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
TBR EXTERNAL HEXAGON Ø3,5	TPDH2L6	2	25	6	2,7	4,5	1,5	2,3	Straight	Semiespherical	1,70 Hexalobular
TBR INTERNAL OCTAGON PLATFORM 3,5	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
TBR INTERNAL OCTAGON PLATFORM 4	TPDH2L74	2	25	7,4	3,3	6	1,4	2,3	Straight	Semiespherical	1,70 Hexalobular
THOMMEN INTERNAL PLATFORM Ø3,5	TPDH14L52	1,4	15	5,2	2,9	3,4	1,8	2,3	Straight	45° Chamfer	1,70 Hexalobular



DYNAMIC CAD-CAM SCREWS DESCRIPTION

IMPLANT	REFERENCE	METRIC	TORQUE N-cm	LENGTH (mm)					HEAD DIAMETER mm	SETTLEMENT	ANGLE	THREAD ENTRY	CONNECTION
				TOTAL	THREAD	A	B	C					
THOMMEN INTERNAL PLATFORM Ø4	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
THOMMEN INTERNAL PLATFORM Ø4,5	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
THOMMEN INTERNAL PLATFORM Ø5	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
THOMMEN INTERNAL PLATFORM Ø6	TPDH16L76	1,6	20	7,6	3,6	6,1		1,5	2,3	Straight		Semiespherical	1,70 Hexalobular
TRI IMPLANTS VENT Ø3,5	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
TRI IMPLANTS TRI-NARROW	TPDH18L80	1,8	25	8	4	6,5		1,5	2,3	Straight		45° Chamfer	1,70 Hexalobular
ZIMMER SCREW-VENT NARROW PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ZIMMER SCREW-VENT REGULAR PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ZIMMER SCREW-VENT WIDE PLATFORM	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ZIMMER SWISS PLUS Ø3,8	TPDH17L71	N1-72	25	7,1	2,6	5,56	5,65	1,45	2,3	Conical	70°	45° Chamfer	1,70 Hexalobular
ZIMMER SWISS PLUS REGULAR PLATFORM	TPDH2L74	2	25	7,4	3,3	6		1,4	2,3	Straight		Semiespherical	1,70 Hexalobular