

ScAnalog
DYNAMIC ABUTMENT® SOLUTIONS



List of compatibilities available*

CODE	MAIN COMPATIBILITY*	REFERENCE	LIBRARY CODE
0008	Astra Evolution Uni Abutment	23.413.008.01-2	DAS_SA_0008
0015	Megagen AnyRidge RP	23.413.015.01-2	DAS_SA_0015
0021	Nobel Biocare Active NP	23.412.021.01-2	DAS_SA_0021
0022	Nobel Biocare Active RP	23.413.022.01-2	DAS_SA_0022
0023	Nobel Biocare Branemark NP	23.412.023.01-2	DAS_SA_0023
0024	Nobel Biocare Branemark RP	23.413.024.01-2	DAS_SA_0024
0025	Nobel Biocare Multi Unit RP	23.413.025.01-2	DAS_SA_0025
0030	Osstem Implant TS RP	23.413.030.01-2	DAS_SA_0030
0040	Zimmer Screw-Vent 3,5	23.412.040.01-2	DAS_SA_0040
0041	Zimmer Screw-Vent 4,5	23.413.041.01-2	DAS_SA_0041
0074	Straumann Synocta 4,8	23.413.074.01-2	DAS_SA_0074
0101	Straumann Screw-Retained NC/RC	23.413.101.01-2	DAS_SA_0101

*Check new compatibilities available emailing:
das@dynamicabutment.com

ScAnalog



AV. BLONDEL, 54 3º
25002 LLEIDA (SPAIN)
TEL. +34 873 450 709
das@dynamicabutment.com
www.dynamicabutment.com



DYNAMIC ABUTMENT® SOLUTIONS

Dynamic Abutment® Solutions is a company with more than 15 years of experience in the research and development of products for the dental sector. It has become a technological expert in angulated solutions and digital treatment of dental prostheses.


The ScAnalog from Dynamic Abutment® Solutions is a scanbody with the dental implant connection that allows the digital recording of the implant position by means of extraoral scan directly on the impression tray.

ScAnalog
DYNAMIC ABUTMENT® SOLUTIONS


3shape

exocad

dental wings



*Case made with ScAnalog compatible with Nobel Biocare Multi Unit RP (Code 0025)



*Case made with ScAnalog compatible with Nobel Biocare Active RP (Code 0022)

