

DOUBLE GINGIVAL HEIGHT


As our main objective is to make your job easier, our R&D&i is committed to meet your needs by developing products that resolve your problems arising from unfavorable placement of implants and improve the aesthetics and functionality of the prosthesis.

For this reason, we have developed a 2mm height Dynamic TiBase® to ease your work in those cases when the gingival margins are at different levels.

Our dynamic TiBase® system provides you with a simple solution for correcting angles up to 30° for both unitary and multiple prosthesis.

Dynamic TiBase® compatible with Straumann Bone Level 1 mm gingival height

Dynamic TiBase® compatible with Straumann Bone Level 2 mm gingival height
Maximum angulation 20°



Angulated
Own Technology
DYNAMIC ABUTMENT SOLUTIONS

Download our new catalogue of CAD-CAM SOLUTIONS 2016 on our website or ask for a free copy to the following e-mail address: das@dynamicabutment.com

 DYNAMIC
ABUTMENT®
SOLUTIONS

AV. BLONDEL, 54 3º · 25002·
LLEIDA (SPAIN) · TEL. +34 973 289 580
das@dynamicabutment.com
www.dynamicabutment.com

