

DYNAMIC TIBASE® 3.0



Compatible with - ASTRA

MODEL	IMPLANT DIAMETER	PLAT.	DYNAMIC TIBASE® 3.0			DYNAMIC SCREW	STRAIGHT SCREW	ANALOG	SCANBODY	LIBRARY CODE
			(R) (NR)	NON-ENGAGING ENGAGING	GINGIVAL HEIGHT					
Evolution 4.8	4,8	4,8	R	31.324.091.01-2	1,2	41.320.082.01-2	40.320.005.01-2	--	30.413.002.01-2	DAS_E_0091
			NR	31.314.091.01-2						
Yellow	3	Yellow (X-narrow)	R	31.322.109.01-2	1,2	41.314.070.01-2	--	--	30.412.001.01-2	DAS_E_0109
			NR	31.312.109.01-2						

Compatible with- BEGO

MODEL	IMPLANT DIAMETER	PLAT.	DYNAMIC TIBASE® 3.0			DYNAMIC SCREW	STRAIGHT SCREW	ANALOG	SCANBODY	LIBRARY CODE
			(R) (NR)	NON-ENGAGING ENGAGING	GINGIVAL HEIGHT					
Semados /SRI 5.5	5,5	5,5	R	31.325.081.01-2	0,4	41.318.064.01-2	--	--	30.414.003.01-2	DAS_E_0081
			NR	31.315.081.01-2						

Compatible with - BIOMET 3i

MODEL	IMPLANT DIAMETER	PLAT.	DYNAMIC TIBASE® 3.0			DYNAMIC SCREW	STRAIGHT SCREW	ANALOG	SCANBODY	LIBRARY CODE
			(R) (NR)	NON-ENGAGING ENGAGING	GINGIVAL HEIGHT					
Certain	5	5	R	31.324.057.01-2	0,5	41.316.084.01-2	40.316.003.01-2	22.614.057.01-2	30.414.003.01-2	DAS_E_0057
			NR	31.314.057.01-2						
Osseotite	5	5	R	31.324.058.01-2	0,5	41.320.047.01-2	40.320.003.01-2	22.614.058.01-2	30.414.003.01-2	DAS_E_0058
			NR	31.314.058.01-2						

Compatible with - BTI

MODEL	IMPLANT DIAMETER	PLAT.	DYNAMIC TIBASE® 3.0			DYNAMIC SCREW	STRAIGHT SCREW	ANALOG	SCANBODY	LIBRARY CODE
			(R) (NR)	NON-ENGAGING ENGAGING	GINGIVAL HEIGHT					
Internal WP	4,5/5/5,5/6/6,5	Wide 5,5	R	31.324.059.01-2	0,7	41.318.065.01-2	40.318.003.01-2	22.614.059.01-2	30.414.003.01-2	DAS_E_0059
			NR	31.314.059.01-2						
External WP	4,5/5/5,5	Wide 5,5	NR	31.314.060.01-2	0,5	41.320.060.01-2	40.320.003.02-2	22.614.060.01-2	30.415.007.01-2	DAS_E_0060

Compatible with - DENTAURUM

MODEL	IMPLANT DIAMETER	PLAT.	DYNAMIC TIBASE® 3.0			DYNAMIC SCREW	STRAIGHT SCREW	ANALOG	SCANBODY	LIBRARY CODE
			(R) (NR)	NON-ENGAGING ENGAGING	GINGIVAL HEIGHT					
Tiologic M	3,7/4,2	Medium	R	31.323.131.01-2	0,5	41.316.081.01-2	--	--	30.413.002.01-2	DAS_E_0131
			NR	31.313.131.01-2						

Compatible with - CAMLOG

MODEL	IMPLANT DIAMETER	PLAT.	DYNAMIC TIBASE® 3.0			DYNAMIC SCREW	STRAIGHT SCREW	ANALOG	SCANBODY	LIBRARY CODE
			(R) (NR)	NON-ENGAGING ENGAGING	GINGIVAL HEIGHT					
Conelog Screw-Line	3,8	3,8	R	31.323.121.01-2	1	41.316.080.01-2	--	--	30.413.002.01-2	DAS_E_0120
			NR	31.313.121.01-2						
	4,3	4,3	R	31.323.121.01-2	1	41.316.080.01-2	--	--	30.413.002.01-2	DAS_E_0121
			NR	31.313.121.01-2						

Compatible with - MEGAGEN

MODEL	IMPLANT DIAMETER	PLAT.	DYNAMIC TIBASE® 3.0			DYNAMIC SCREW	STRAIGHT SCREW	ANALOG	SCANBODY	LIBRARY CODE
			(R) (NR)	NON-ENGAGING ENGAGING	GINGIVAL HEIGHT					
Octógono Externo 3.8	Converter AnyRidge/AnyOne	3,8	R	31.322.128.01-2	2,5	41.320.044.01-2	40.320.002.01-2	--	30.413.002.01-2	DAS_E_0128

Compatible with - NOBEL BIOCARE

MODEL	IMPLANT DIAMETER	PLAT.	DYNAMIC TIBASE® 3.0			DYNAMIC SCREW	STRAIGHT SCREW	ANALOG	SCANBODY	LIBRARY CODE
			(R) (NR)	NON-ENGAGING ENGAGING	GINGIVAL HEIGHT					
Active	5,5	Wide	R	31.324.124.01-2	1,4	41.320.075.01-2	40.320.008.02-2	--	30.414.003.01-2	DAS_E_0124
			NR	31.314.124.02-2						
Branemark WP	5/6	Wide	R	31.324.061.01-02	0,6	41.325.067.01-2	--	22.614.061.01-2	30.415.007.01-2	DAS_E_0061
			NR	31.314.061.01-02						