



DYNAMIC ABUTMENT • SOLUTIONS

# MULTI-UNIT DASSYSTEM

Work different, work better

## LIBRARY CODES

DYNAMIC ABUTMENT® SOLUTIONS LIBRARY



**NOMENCLATURE:**

DAS = Dynamic Abutment Solutions    SA = ScAnalog  
 MU = Multi-Unit    HH = Scanbody height, can be 8mm, 10mm or 12mm  
 XXXX = Compatibility (Example 0040)    I/IG = Standard Adapter I or Adapter 3mm IG  
 PREMIUM = It has the prosthetic options of Tibases, Direct to Implant, Direct to Provisional Impl both Direct to Implant with real connection values for milling. Also the option of Digital Analog.

4,5 = MU Scanbody height (mm)  
 8 = MU Scanbody height (mm)

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
BIOMET 3I	Certain	3,25/4	3,4	DAS_MU_0001
BIOMET 3I	Certain	4/5	4,1	DAS_MU_0002
AB	I2	3,5/3,75/4,2/4,5/5/6	Standard	DAS_MU_0040
AB	I22	3,75/4,2	Standard	DAS_MU_0040
AB	I5	3,5/3,75/4,2/4,5/5/6/7/8	Standard	DAS_MU_0040
AB	I55	3,75/4,2/4,5/5/6/7/8	Standard	DAS_MU_0040
AB	I10	4,2/5	Standard	DAS_MU_0040
AB	I15	6/7/8	Standard	DAS_MU_0040
ACE	Infinity Internal Hex	3,7/4,1	3,5	DAS_MU_0040
ADIN	Swell	3,3	3,45	DAS_MU_0040
ADIN	Swell	3,75/4,2	3,6	DAS_MU_0040
ADIN	Swell	5	4	DAS_MU_0040
ADIN	Swell	6	4,6	DAS_MU_0040
ADIN	Toureg-CloseFit	3,5	RP	DAS_MU_0021
ADIN	Toureg-CloseFit	4,3/5	WP	DAS_MU_0022
ADIN	Toureg-S/Toureg-OS	3,5	3,45	DAS_MU_0040
ADIN	Toureg-S/Toureg-OS	3,75/4,2	3,6	DAS_MU_0040
ADIN	Toureg-S/Toureg-OS	5	4	DAS_MU_0040
ADIN	Toureg-S/Toureg-OS	6	5	DAS_MU_0040
ALFA-GATE	Bioactive S-line/SLA	3,3/3,75/4,2/4,7/5,2/6	3,75	DAS_MU_0040
ALPHABIO	Internal Hex Connection (IH) SPI	3,3/3,75/4,2/5/6	Universal	DAS_MU_0040
ALPHABIO	Internal Hex Connection (IH) ICE	3,7/3,75/4,2/4,65/5,3	Universal	DAS_MU_0040
ALPHABIO	Internal Hex Connection (IH) DFI	3,3/3,75/4,2/4,5	Universal	DAS_MU_0040
ALPHABIO	Internal Hex Connection (IH) NEO	3,75/4,2/5	3,5	DAS_MU_0040
ALPHABIO	Conical Standard Connection (CS)	3,75/4,2/5	Standard	DAS_MU_0169
ANCLADEN	Anclalock	3,75/4,25/5	3,5	DAS_MU_0040
ANKYLOS	Ankylos		3,5/4,5/5,5/7	DAS_MU_0075
ARDS	Classic	3,75/4,2/4,5	3,75	DAS_MU_0040
ARDS	Classic	3,3/3,75/4,2/5/6	3,75	DAS_MU_0040
ARDS	CIT	3,3/3,75/4,2/5/6	3,75	DAS_MU_0040
ARDS	Premium	3,3/3,75/4,2/5/6	3,75	DAS_MU_0040

**NOMENCLATURE:**

DAS = Dynamic Abutment Solutions    SA = ScAnalog    4,5 = MU Scanbody height (mm)  
 MU = Multi-Unit    HH = Scanbody height, can be 8mm, 10mm or 12mm    height (mm)  
 XXXX = Compatibility (Example 0040)    I/IG = Standard Adapter I or Adapter 3mm IG    8 = MU Scanbody height (mm)  
 PREMIUM = It has the prosthetic options of Tibases, Direct to Implant, Direct to Provisional Implant, both Direct to Implant with real connection values for milling. Also the option of Digital Analog.



BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
ARDS	Smart	3,75/4,2/4,5	3,75	DAS_MU_0040
ASTRA	Aqua	3,5/4	Aqua(Estrecha)	DAS_MU_0004
ASTRA	Evolution	3,6	3,6	DAS_MU_0006
ASTRA	Evolution	4,2	4,2	DAS_MU_0007
ASTRA	Evolution	4,8	4,8	DAS_MU_0091
AVINENT	HI	3,1//3,5/4	3,5	DAS_MU_0040
AVINENT	HI	3,3/3,8/4,4/2,4,8//4,5/5	4,1	DAS_MU_0040
BIOCONCEPT	BC Bone Level	3,3	Narrow	DAS_MU_0033
BIOCONCEPT	BC Bone Level	4,1/4,8	Regular	DAS_MU_0035
BIOCONCEPT	BV Tapered Bone Level	3,5	Narrow	DAS_MU_0029
BIOCONCEPT	BV Tapered Bone Level	4/4,5/5	Regular	DAS_MU_0030
BIOHORIZON	Internal	3/3,8/4,6	3,5 (Yellow)	DAS_MU_0040
BIOHORIZON	Tapered Tissue Level	3/3,8	3,5 (Yellow)	DAS_MU_0040
BIOHORIZON	Hexágono Interno	3/3,4	3	DAS_MU_0102
BIONER	Hikelt	3,8	3,95	DAS_MU_0040
BIONER	TopDM	3,5	3,5	DAS_MU_0021
BIONER	TopDM	4	4	DAS_MU_0021
BIONER	TopDM	5	5	DAS_MU_0021
BIOTEC	SPR/CIM	3,3	3,3	DAS_MU_0040
BIOTEC	SPR/CIM	3,75	3,75	DAS_MU_0040
BIOTEC	SPR/SPTT/CIM	4,2	4,2	DAS_MU_0040
BIOTEC	SPR/SPTT/CIM	5	5	DAS_MU_0040
BTK	Klassic / Konic INT	3,25/4	3,5 IR	DAS_MU_0040
B&W	Cónico Hexagono Interno CIH	3,3/4	4	DAS_MU_0040
CAMLOG	Conelog Screw-Line	3,8	3,8	DAS_MU_0120
CAMLOG	Conelog Screw-Line	4,3	4,3	DAS_MU_0121
CORTEX	Classix	3,3/3,8/4,2/5/6	3,75	DAS_MU_0040
CORTEX	Dynamix	3,3/3,8/4,2/5/6	3,75	DAS_MU_0040
CORTEX	Saturn	3,8/4,2	3,5	DAS_MU_0040
CORTEX	Conical Platform	3,3/3,8/4,2	RP	DAS_MU_0004
CORTEX	Magix	3,3/3,8/4,2	RP	DAS_MU_0004

**NOMENCLATURE:**

DAS = Dynamic Abutment Solutions    SA = ScAnalog  
 MU = Multi-Unit    HH = Scanbody height, can be 8mm, 10mm or 12mm  
 XXXX = Compatibility (Example 0040)    I/IG = Standard Adapter I or Adapter 3mm IG  
 PREMIUM = It has the prosthetic options of Tibases, Direct to Implant, Direct to Provisional Impl both Direct to Implant with real connection values for milling. Also the option of Digital Analog.

4,5 = MU Scanbody height (mm)  
 8 = MU Scanbody height (mm)

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
CONEXÃO SISTEMA DE PRÓTESE	Flash	3,5/4,3/5	Universal	DAS_MU_0021
CONEXÃO SISTEMA DE PRÓTESE	Torq	3,5/3,75/4	Universal	DAS_MU_0021
CONEXÃO SISTEMA DE PRÓTESE	Expand	3,75/4/5	Universal	DAS_MU_0021
DENTIUM	SuperLine-Implantium	3,4	3,6	DAS_MU_0030
DENTIUM	SuperLine-Implantium	3,8	4	DAS_MU_0030
DENTIUM	SuperLine-Implantium	4,3	4,5	DAS_MU_0030
DENTIUM	SuperLine-Implantium	4,8	5	DAS_MU_0030
DENTIUM	SuperLine-Implantium	4,8	6	DAS_MU_0030
DENTIS	S-Clean	3,7	Mini	DAS_MU_0030
DENTIS	S-Clean	4,1/4,3	Regular	DAS_MU_0030
DENTIS	S-Clean	4,8	Wide	DAS_MU_0030
DIO IMPLANTS	UF	3,8/4/4,5/5/5,5	Regular	DAS_MU_0030
EASY IMPLANT	Master C	3,5	3,5 (Ocean)	DAS_MU_0004
EASY IMPLANT	Master C	4	4 (Ocean)	DAS_MU_0004
EASY IMPLANT	Master C	4,5	4,5 (Lilas)	DAS_MU_0030
EASY IMPLANT	Master C	5	5 (Lilas)	DAS_MU_0030
EASY IMPLANT	Master S	3,3	3,3 (Ocean)	DAS_MU_0004
EASY IMPLANT	Master S	3,75	3,75 (Lilas)	DAS_MU_0030
EASY IMPLANT	Master S	4,25	4,25 (Lilas)	DAS_MU_0030
EASY IMPLANT	Master S	4,75	4,75 (Lilas)	DAS_MU_0030
EASY IMPLANT	Master L	3,3	3,3 (Lilas)	DAS_MU_0030
EASY IMPLANT	Master L	3,75	3,75 (Lilas)	DAS_MU_0030
EASY IMPLANT	Master L	4,25	4,25 (Lilas)	DAS_MU_0030
EASY IMPLANT	Master L	4,75	4,75 (Lilas)	DAS_MU_0030
EUROTEKNIKA	Naturactis	3,5	3,4	DAS_MU_0004
EUROTEKNIKA	Naturactis	4	3,8	DAS_MU_0004
EUROTEKNIKA	Naturactis	4,5	4,3	DAS_MU_0004
EUROTEKNIKA	Naturactis	5	4,8	DAS_MU_0004
EUROTEKNIKA	Natea	3,6/4,1/4,8	Narrow	DAS_MU_0004
EUROTEKNIKA	Natea	3,6/4,1/4,8	Regular	DAS_MU_0004
EUROTEKNIKA	Naturall	3,5	Narrow	DAS_MU_0004

**NOMENCLATURE:**

DAS = Dynamic Abutment Solutions    SA = ScAnalog  
 MU = Multi-Unit    HH = Scanbody height, can be 8mm, 10mm or 12mm  
 XXXX = Compatibility (Example 0040)    I/IG = Standard Adapter I or Adapter 3mm IG  
 PREMIUM = It has the prosthetic options of Tibases, Direct to Implant, Direct to Provisional Implant, both Direct to Implant with real connection values for milling. Also the option of Digital Analog.

4,5 = MU Scanbody height (mm)  
 8 = MU Scanbody height (mm)

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
EUROTEKNIKA	Naturall	4/4,5	Regular	DAS_MU_0004
GALIMPLANT	Conexióninterna	3,5	3,5	DAS_MU_0004
GALIMPLANT	Conexióninterna	4	4	DAS_MU_0004
GALIMPLANT	Conexióninterna	5	5	DAS_MU_0004
GMI	Frontier	3,3/3,75/4,25	RP 3,3	DAS_MU_0040
HAHN IMPLANT	Tapered Implant	3,5/4,3	Purple	DAS_MU_0021
HAHN IMPLANT	Hahn Tapered Implant	5	5	DAS_MU_0022
HIOSSSEN	ETII SA / ETIII SA	3,5	Mini	DAS_MU_0029
HIOSSSEN	ETII SA / ETIII SA	4/4,5/5	Regular	DAS_MU_0030
HIOSSSEN	ETIII BA	3,5	Mini	DAS_MU_0029
HIOSSSEN	ETIII BA	4/4,5/5	Regular	DAS_MU_0030
HI-TEC	Tapered Self Thread	3,3/3,75	3,5	DAS_MU_0040
HI-TEC	Logic Plus	3,5	3,7	DAS_MU_0040
HI-TEC	Logic Plus	4,3	3,9	DAS_MU_0040
IBS	Magic FC	4/4,5/5/5,5/6/6,5	3,8	DAS_MU_0030
IBS	N.R. Fix	3/3,5	3,8	DAS_MU_0030
IDO IMPLANTS	I Do Implant	3,8/4/4,5/5/5,5/6/7	Universal	DAS_MU_0030
IMPLANT DIRECT	Interactive	3,2/3,7	3	DAS_MU_0021
IMPLANT DIRECT	Interactive	4,3/5	3,4	DAS_MU_0022
IMPLANT DIRECT	Legacy	3,7	3,5	DAS_MU_0040
IMPLANT DIRECT	SwishActive	3,3	3	DAS_MU_0021
IMPLANT DIRECT	SwishActive	4,1/4,8	3,4	DAS_MU_0022
IMPLANT GENESIS	Aktive System	3,5/3,75/4,2/5	Standard	DAS_MU_0040
IHDE DENTAL (IMBIODENT)	Bone Level Plus	3,3	3,3	DAS_MU_0033
IHDE DENTAL (IMBIODENT)	Bone Level Plus	4,1	4,1	DAS_MU_0035
IHDE DENTAL (IMBIODENT)	Bone Level Plus	4,1	4,8	DAS_MU_0035
INTRA LOCK	IntraHex	3,75/4	3,5	DAS_MU_0040
JDENTALCARE	JDEvolution/JDEvolution Plus	3,7	3,7	DAS_MU_0040
JDENTALCARE	JDEvolution/JDEvolution Plus	4,3/5	4	DAS_MU_0040
JDENTALCARE	JDEvolution/JDEvolution Plus	6	5	DAS_MU_0040
JDENTALCARE	JD ICON	3,9	3,9	DAS_MU_0022

**NOMENCLATURE:**

DAS = Dynamic Abutment Solutions

SA = ScAnalog

MU = Multi-Unit

HH = Scanbody height, can be 8mm, 10mm or 12mm

XXXX = Compatibility (Example 0040)

I/IG = Standard Adapter I or Adapter 3mm IG

PREMIUM = It has the prosthetic options of Tibases, Direct to Implant, Direct to Provisional Impl both Direct to Implant with real connection values for milling. Also the option of Digital Analog.

4,5 = MU Scanbody height (mm)

8 = MU Scanbody height (mm)

**NOMENCLATURE:**

DAS = Dynamic Abutment Solutions

SA = ScAnalog

MU = Multi-Unit

HH = Scanbody height, can be 8mm, 10mm or 12mm

XXXX = Compatibility (Example 0040)

I/IG = Standard Adapter I or Adapter 3mm IG

PREMIUM = It has the prosthetic options of Tibases, Direct to Implant, Direct to Provisional Implant, both Direct to Implant with real connection values for milling. Also the option of Digital Analog.

4,5 = MU Scanbody height (mm)

8 = MU Scanbody height (mm)



BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
JDENTALCARE	JD ICON	4,3	4	DAS_MU_0022
JDENTALCARE	JD ICON	5	4,7	DAS_MU_0022
LEADER	Tixos Internal Hex	3,3	3,5	DAS_MU_0040
LEADER	Tixos Internal Hex	3,75	4	DAS_MU_0040
MEDENTIS	ICX-Templant	3,75	3,75	DAS_MU_0125
MEDENTIS	ICX-Templant	4,1	4,1	DAS_MU_0125
MEDENTIS	ICX-Templant	4,8	4,8	DAS_MU_0125
MEGAGEN	AnyOne Internal	3,5/4/4,5/5/6/7	General	DAS_MU_0030
MEGAGEN	AnyRidge	3,5	Small	DAS_MU_0015
MEGAGEN	AnyRidge	4/4,5	Regular	DAS_MU_0015
MEGAGEN	AnyRidge	5/5,5	Wide	DAS_MU_0015
MICRODENT	Ektos	3,7/4,2	3,5	DAS_MU_0040
MIS	M4	3,75/4,2	Standard	DAS_MU_0040
MIS	Seven	3,75/4,2	Standard	DAS_MU_0040
MOZO-GRAU	MG Inhex	3,75/4,25	2,8 Standard	DAS_MU_0004
MPI	Excellence CM	3,5/4	Regular	DAS_MU_0004
NEOBIOTECH	IS System	4	Regular 4	DAS_MU_0030
NEOBIOTECH	IS System	4,5	Regular 4,5	DAS_MU_0030
NEOBIOTECH	IS System	5	Wide 5	DAS_MU_0030
NEODENT	GM			DAS_MU_0186
NOBEL BIOCARE	Active	3,5	Narrow	DAS_MU_0021
NOBEL BIOCARE	Active	4,3/5	Regular	DAS_MU_0022
NORIS MEDICAL	TUFF	3,3/3,75/4,2/5/6	3,75	DAS_MU_0040
NORIS MEDICAL	TUFF TT	3,3/3,75/4,2/5/6	3,75	DAS_MU_0040
NORIS MEDICAL	ONIX	3,3/3,75/4,2/5/6	3,75	DAS_MU_0040
NORIS MEDICAL	CORTICAL	4/5/6.	3,75	DAS_MU_0040
NORIS MEDICAL	PTERYCORE	4,2	3,75	DAS_MU_0040
NORIS MEDICAL	PTERYFIT	4,2	3,75	DAS_MU_0040
NORMON	HI	3,75/4,25/4,75	3,5	DAS_MU_0040
NOVA IMPLANTS	PSI/PCI	3,3/3,75/4,2/5/6	3,75	DAS_MU_0040
OSSTEM	TS	3,5	Mini	DAS_MU_0029

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
OSSTEM	TS	4/4,5/5/6/7	Regular	DAS_MU_0030
OSTEOPLUS	SHI	3,3 / 3,75 / 4,2	3,5	DAS_MU_0040
PALTOP	Conical Active	3,75/4,2/5	Standard	DAS_MU_0029
OXY	K1 Line	3,5/4/4,5/5	Regular	DAS_MU_0015
OXY	PSK Line	3,5/4/4,5/5	Regular	DAS_MU_0015
OXY	MD Line KONE	4,75/4,25/5	Regular	DAS_MU_0015
PROCLINIC	Cilindrico Interno/Cónico Interno		3,5	DAS_MU_0040
PROCLINIC	Aqua CM	3,5/4/5	2,82	DAS_MU_0004
SEWON MEDIX	IH2 SLA SYSTEM	3,5	Mini	DAS_MU_0029
SEWON MEDIX	IH2 RBM SYSTEM	3,5	Mini	DAS_MU_0029
SEWON MEDIX	IH2 SLA SYSTEM	3,5/4/4,5/5	Regular	DAS_MU_0030
SEWON MEDIX	IH2 RBM SYSTEM	3,5/4/4,5/5	Regular	DAS_MU_0030
SGS	P1	3,2/3,75/4,2/5/6	3,75	DAS_MU_0040
SGS	P7	3,2/3,75/4,2/4,5/5/6	3,75	DAS_MU_0040
SIGNO VINCES	Compact	4,5	CM3,8	DAS_MU_0004
SIGNO VINCES	Duocon	3,8	CM3,8	DAS_MU_0004
SIGNO VINCES	Infra	3,3/3,8/4,6	CM	DAS_MU_0004
SOUTHERN IMPLANTS	Internal Hex	3,75/4,2/5	Universal	DAS_MU_0040
SOUTHERN IMPLANTS	Deep Conical	3,5/4	2,95/3,1	DAS_MU_0004
STRAUMANN	Bone Level	3,3	NC- 3,3	DAS_MU_0033
STRAUMANN	Bone Level	4,1	RC-4,1	DAS_MU_0035
STRAUMANN	BLX	3,5/3,75/4/4,5	RB	DAS_MU_0207
TITANIUM-FIX	B-fix	3,5/4	Regular	DAS_MU_0004
TRI DENTAL IMPLANTS	Tri-Vent	3,75/4,1/4,7	3,5	DAS_MU_0040
TREE-OSS	Simple	3,3/3,75/5	3,75	DAS_MU_0040
UFIT	Gt2	3,5	Mini	DAS_MU_0004
UFIT	Nt2	3,5	Mini	DAS_MU_0004
YES IMPLANTS	S-SYSTEM	3/3,5	Narrow	DAS_MU_0030
YES IMPLANTS	S-SYSTEM	4/4,5	Regular	DAS_MU_0030
YES IMPLANTS	S-SYSTEM	5/5,5	Wide	DAS_MU_0030
YES IMPLANTS	S-SYSTEM	6/6,5/7/8	B-Wide	DAS_MU_0030

**NOMENCLATURE:**

DAS = Dynamic Abutment Solutions

SA = ScAnalog

MU = Multi-Unit

HH = Scanbody height, can be 8mm, 10mm or 12r

XXXX = Compatibiliy (Example 0040)

I/IG = Standard Adapter I or Adapter 3mm IG

**PREMIUM** = It has the prosthetic options of Tibases, Direct to Implant, Direct to Provisional Impl both Direct to Implant with real connection values for milling. Also the option of Digital Analog.

4,5 = MU Scanbody height (mm)

8 = MU Scanbody height (mm)

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
ZIMMER	Zimmer Type (Internal Hexagon)	3,7/4,1	3,5	DAS_MU_0040
ZIMMER	Screw-Vent	3,7/4,1	3,5	DAS_MU_0040