




**DYNAMIC
ABUTMENT®
SOLUTIONS**

**FULL SYSTEM
(INTRAORAL AND LAB)
LIBRARY CODES**

DYNAMIC ABUTMENT® SOLUTIONS LIBRARY



NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging



BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
BIOMET 3I	Certain	3,25/4	3,4	DAS_I_0001
BIOMET 3I	Certain	4/5	4,1	DAS_I_0002
BIOMET 3I	Certain	5	5	DAS_I_0057
BIOMET 3I	External Hex	3,25	3,4	DAS_I_0003
BIOMET 3I	External Hex	3,75/4	4,1	DAS_I_0024
BIOMET 3I	External Hex	5	5	DAS_I_0058
BIOMET 3I	Low profile	Universal	Universal	DAS_I_0025
AB	I2	3,5/3,75/4,2/4,5/5/6	Standard	DAS_I_0040
AB	I22	3,75/4,2	Standard	DAS_I_0040
AB	I5	3,5/3,75/4,2/4,5/5/6/7/8	Standard	DAS_I_0040
AB	I55	3,75/4,2/4,5/5/6/7/8	Standard	DAS_I_0040
AB	I10	4,2/5	Standard	DAS_I_0040
AB	I15	6/7/8	Standard	DAS_I_0040
AB	Multi Unit D1P64	Universal	Universal	DAS_I_0025
ACE	External Hex	3,3	NP 3,5	DAS_I_0023
ACE	External Hex	3,75/4	RP 4,1	DAS_I_0024
ACE	External Hex	4,75	WP 5	DAS_I_0058
ACE	Infinity Internal Hex	3,7/4,1	3,5	DAS_I_0040
ACE	Infinity Internal Hex	4,7/5,1	4,5	DAS_I_0041
ACE	Infinity TRICAM	3,5	3,5	DAS_I_0026
ACE	Infinity TRICAM	4,3	4,3	DAS_I_0027
ACE	Infinity TRICAM	5	5	DAS_I_0028
ACE	Infinity Octagon	3,3/4,1/4,8	RP 4,8	DAS_I_0037
ACE	Infinity Octagon	4,8	WP 6,5	DAS_I_0096
ACE	Multi Unit	Universal	Universal	DAS_I_0025
ADIN	Swell	3,3	3,45	DAS_I_0040
ADIN	Swell	3,75/4,2	3,6	DAS_I_0040
ADIN	Swell	5	4	DAS_I_0040
ADIN	Swell	6	4,6	DAS_I_0040
ADIN	TouregCloseFit	2,75	UNP	DAS_I_0188
ADIN	TouregCloseFit	3	NP	DAS_I_0145
ADIN	TouregCloseFit	3,5	RP	DAS_I_0021
ADIN	TouregCloseFit	4,3/5	WP	DAS_I_0022

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
ADIN	TouregS/TouregOS	3,5	3,45	DAS_I_0040
ADIN	TouregS/TouregOS	3,75/4,2	3,6	DAS_I_0040
ADIN	TouregS/TouregOS	5	4	DAS_I_0040
ADIN	TouregS/TouregOS	6	5	DAS_I_0040
ADIN	Multi Unit TMA	Universal	Universal	DAS_I_0025
ALFA-GATE	Bioactive S-line/SLA	3,3/3,75/4,2/4,7/5,2/6	3,75	DAS_I_0040
ALFA-GATE	Conical	3,5	NP	DAS_I_0021
ALFA-GATE	Conical	4,3/5	RP	DAS_I_0022
ALPHABIO	Conical Hex Connection (NICE)	3,2	Narrow	DAS_I_0136
ALPHABIO	Conical Hex Connection (NEO)	3,2/3,5	Narrow	DAS_I_0136
ALPHABIO	Internal Hex Connection (IH) SPI	3,3/3,75/4,2/5/6	Universal	DAS_I_0040
ALPHABIO	Internal Hex Connection (IH) ICE	3,7/3,75/4,2/4,65/5,3	Universal	DAS_I_0040
ALPHABIO	Internal Hex Connection (IH) DFI	3,3/3,75/4,2/5	Universal	DAS_I_0040
ALPHABIO	Internal Hex Connection (IH) ATID	3,3/3,75/4,2/5/6	Universal	DAS_I_0040
ALPHABIO	Internal Hex Connection (IH) NEO	3,75/4,2/5	3,5	DAS_I_0040
ALPHABIO	Conical Standard Connection (CS)	3,75/4,2/5	Standard	DAS_I_0169
ANCLADEN	Anclalock	3,75/4,25/5	3,5	DAS_I_0040
ANKYLOS	Ankylos	3,5/4,5/5,5/7	NP/RP/WP	DAS_I_0075
ANKYLOS	Balance Base		Narrow	DAS_I_0183
ANTHOGRYR	Axiom REG / PX	3,4	3,4 XNP	DAS_I_0161
ANTHOGRYR	Axiom REG / PX	4	4 RP	DAS_I_0149
ANTHOGRYR	Axiom REG / PX	4,6	4,6 WP	DAS_I_0149
ANTHOGRYR	Axiom REG / PX	5,2	5,2	DAS_I_0162
ANTHOGRYR	Anthofit HE	3,5/3,75/4	R (4,1)	DAS_I_0024
ANTHOGRYR	Ossfit	3,5/4,2	4,8	DAS_I_0037
ANTHOGRYR	Ossfit	3,5/4,2	4,8	DAS_I_0074
ANTHOGRYR	Ossfit	5	6,5	DAS_I_0096
ANTHOGRYR	Infinity Octagon	4,8	WP 6,5	DAS_I_0096
ANTHOGRYR	Multi Unit	4,8	Universal	DAS_I_0163
ANTHOGRYR	Axiom REG / PX	5,2	5,2	DAS_I_0162
ARDS	Classic	3,75/4,2/4,5	3,75	DAS_I_0040
ARDS	Classic	3,3/3,75/4,2/5/6	3,75	DAS_I_0040
ARDS	CIT	3,3/3,75/4,2/5/6	3,75	DAS_I_0040

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging



BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
ARDS	Premium	3,3/3,75/4,2/5/6	3,75	DAS_I_0040
ARDS	Smart	3,75/4,2/4,5	3,75	DAS_I_0040
ASTRA	Yellow	3	Yellow	DAS_I_0109
ASTRA	Aqua	3,5/4	Aqua(Estrecha)	DAS_I_0004
ASTRA	Lilac	4,5/5	Lilac (Ancha)	DAS_I_0005
ASTRA	Evolution	3	3	DAS_I_0090
ASTRA	Evolution	3,6	3.6	DAS_I_0006
ASTRA	Evolution	4,2	4,2	DAS_I_0007
ASTRA	Evolution	4,8	4,8	DAS_I_0091
ASTRA	Evolution	5,4	5,4	DAS_I_0092
ASTRA	Evolution	Universal	Universal	DAS_I_0008
AVINENT	HE	3,3//3,5/4	3,5	DAS_I_0023
AVINENT	HE	3,3/3,8/4/4,2/4,8//4,5/5	4,1	DAS_I_0024
AVINENT	HI	3,1//3,5/4	3,5	DAS_I_0040
AVINENT	HI	3,3/3,8/4/4,2/4,8//4,5/5	4,1	DAS_I_0040
AVINENT	Transepitelial	Universal	Regular Platform	DAS_I_0025
AVINENT	HE/EC	4,8	5,1	DAS_I_0061
BEGO	RS/RSX	3	3	DAS_I_0049
BEGO	S/RI/RS/RSX	3,25/3,75	3,25/3,75	DAS_I_0050
BEGO	S/RI/RS/RSX	4,1	4,1	DAS_I_0051
BEGO	S/RI/RS/RSX	4,5	4,5	DAS_I_0052
BEGO	S/RI/RS/RSX	5,5	5,5	DAS_I_0081
BEGO	Mini	2,7/2,9/3,1	Mini	DAS_I_0187
BEGO	MultiPlus	Universal	Universal	DAS_I_0150
BIOCONCEPT	BC Tissue Level Standard	3,3/4,1/4,8	Regular	DAS_I_0037
BIOCONCEPT	BC Tissue Level Standard Plus	4,8	Regular	DAS_I_0037
BIOCONCEPT	BC Tissue Level Tapered Effect	4,8	Regular	DAS_I_0037
BIOCONCEPT	BC Bone Level	3,3	Narrow	DAS_I_0033
BIOCONCEPT	BC Bone Level	4,1/4,8	Regular	DAS_I_0035
BIOCONCEPT	BV Tapered Bone Level	3,5	Narrow	DAS_I_0029
BIOCONCEPT	BV Tapered Bone Level	4/4,5/5	Regular	DAS_I_0030
BIOGENESIS	3icon	3,3	Mini (Pink)	DAS_I_0023

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
BIOGENESIS	3icon	3,75/4/4,3/4,5	Regular (Blue)	DAS_I_0024
BIOGENESIS	Aticon	3,5//4/4,5/5	Pink/Blue	DAS_I_0005
BIOGENESIS	Iticon	3,5/4,1/4,8	4,8	DAS_I_0037
BIOHORIZON	External	3,5	3,7 (Yellow)	DAS_I_0023
BIOHORIZON	Tapered Internal	3/3,4	3.0 (Grey)	DAS_I_0102
BIOHORIZON	Tapered Internal	3,8	3,5 (Yellow)	DAS_I_0040
BIOHORIZON	Tapered Internal	4,6	4,5 (Green)	DAS_I_0041
BIOHORIZON	Tapered Internal	5,8	5,7 (Blue)	DAS_I_0080
BIOHORIZON	Tapered Plus	3,8	3.0 (Grey)	DAS_I_0102
BIOHORIZON	Tapered Plus	4,6	3,5 (Yellow)	DAS_I_0040
BIOHORIZON	Tapered Plus	5,8	4,5 (Green)	DAS_I_0041
BIOHORIZON	Mount-free Tapered Internal	3/3,4	3.0 (Grey)	DAS_I_0102
BIOHORIZON	Mount-free Tapered Internal	3,8	3,5 (Yellow)	DAS_I_0040
BIOHORIZON	Mount-free Tapered Internal	4,6	4,5 (Green)	DAS_I_0041
BIOHORIZON	Mount-free Tapered Internal	5,8	5,7 (Blue)	DAS_I_0080
BIOHORIZON	Tapered Tissue Level	3/3,8	3,5 (Yellow)	DAS_I_0040
BIOHORIZON	Tapered Tissue Level	4,6	4,5 (Green)	DAS_I_0041
BIOHORIZON	Tapered Tissue Level	5,8	5,7 (Blue)	DAS_I_0080
BIOHORIZON	Tapered Pro	3,8	3.0 (Grey)	DAS_I_0102
BIOHORIZON	Tapered Pro	4,2/4,6	3,5 (Yellow)	DAS_I_0040
BIOHORIZON	Tapered Pro	5,2	4,5 (Green)	DAS_I_0041
BIOHORIZON	Tapered Short	4,6	3,5 (Yellow)	DAS_I_0040
BIOHORIZON	Tapered Short	5,8	4,5 (Green)	DAS_I_0041
BIOHORIZON	Tapered TPG	4,2	3,5 (Yellow)	DAS_I_0040
BIOHORIZON	Tapered IM (Inmediate Molar)	7/8	5,7 (Blue)	DAS_I_0080
BIOHORIZON	Multi Unit Angled	Universal	Universal	DAS_I_0025
BIOLOK	External Hex	3,45	3,45	DAS_I_0003
BIONER	Hikelt	3,8	3,95	DAS_I_0040
BIONER	Hikelt	4,7	4,9	DAS_I_0041
BIONER	Ikelt/Bikelt	3,3/3,75/4	4,1	DAS_I_0024
BIONER	Transepitelial A5M	Universal	Regular	DAS_I_0025
BIONER	TopDM	3,5	3,5	DAS_I_0021
BIONER	TopDM	4	4	DAS_I_0021

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging



BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
BIONER	TopDM	5	5	DAS_I_0021
BIOTEC	SPR/CIM	3,3	3,3	DAS_I_0040
BIOTEC	SPR/CIM	3,75	3,75	DAS_I_0040
BIOTEC	SPR/SPTT/CIM	4,2	4,2	DAS_I_0040
BIOTEC	SPR/SPTT/CIM	5	5	DAS_I_0040
BIOTECH	Kontakt	3	Universal	DAS_I_0164
BIOTECH	Kontakt	3,6/4,2/4,8/5,4	Universal	DAS_I_0165
BREDDENT MEDICAL	Narrow Sky	3,5	NP 3,5	DAS_I_0110
BREDDENT MEDICAL	Blue Sky	3,5/4/4,5/5,5	4	DAS_I_0111
BREDDENT MEDICAL	Blue Sky Classic	3,5/4/4,5	4	DAS_I_0111
BTI	Externa Tiny	2,5/3/3,3/3,5/3,75	Tiny 3,5	DAS_I_0009
BTI	Externa Universal	3,75/4	Universal 4,1	DAS_I_0024
BTI	Externa Universal Plus	4,5/5	Universal Plus 4,1	DAS_I_0024
BTI	External Connection	4,5/5/5,5	Ancha 5,5	DAS_I_0060
BTI	Interna Universal	3,3/3,5/3,75/4/4,25	Universal 4,1	DAS_I_0010
BTI	Interna Universal Plus	4,5/5/5,5	Universal Plus 4,1	DAS_I_0010
BTI	Interna ancha	5,5/6/6,25	Ancha 5,5	DAS_I_0059
BTI	Multi IM	Universal	Universal	DAS_I_0151
BTK	Klassic / Konic EXT	3,25	3,4 EN	DAS_I_0003
BTK	Klassic / Konic EXT	3,25PL/3,75/4	4,1 ER	DAS_I_0024
BTK	Klassic / Konic INT	3,25/3,75/4	3,5 IR	DAS_I_0040
B&W	Internal Connection	3,3/4	4	DAS_I_0040
B&W	External Connection	3,75/4	4,1	DAS_I_0024
CAMLOG	Camlog Progressive-Line	3,3	3,3	DAS_I_0087
CAMLOG	Camlog Progressive-Line	3,8	3,8	DAS_I_0011
CAMLOG	Camlog Progressive-Line	4,3	4,3	DAS_I_0012
CAMLOG	Camlog Progressive-Line	5	5	DAS_I_0088
CAMLOG	Camlog Screw-Line	3,3	3,3	DAS_I_0087
CAMLOG	Camlog Screw-Line	3,8	3,8	DAS_I_0011
CAMLOG	Camlog Screw-Line	4,3	4,3	DAS_I_0012
CAMLOG	Camlog Screw-Line	5	5	DAS_I_0088
CAMLOG	Conelog Progressive-Line	3,3	3,3	DAS_I_0119
CAMLOG	Conelog Progressive-Line	3,8	3,8	DAS_I_0120

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
CAMLOG	Conelog Progressive-Line	4,3	4,3	DAS_I_0121
CAMLOG	Conelog Progressive-Line	3,3	3,3	DAS_I_0119
CAMLOG	Conelog Screw-Line	3,8	3,8	DAS_I_0120
CAMLOG	Conelog Screw-Line	4,3	4,3	DAS_I_0121
CORTEX	Conical Dynamix	3	NP	DAS_I_0109
CORTEX	Conical Dynamix	3,3/3,5/4,2	RP	DAS_I_0004
CORTEX	Conical Dynamix	5/6	WP	DAS_I_0005
CORTEX	Conical Classix	3,3/3,8/4,2	RP	DAS_I_0004
CORTEX	Conical Classix	5/6	WP	DAS_I_0005
CORTEX	Conical Magix	3,3/3,8/4,2	RP	DAS_I_0004
CORTEX	Internal Hex Dynamix	3,3/3,8/4,2/5/6	3,5	DAS_I_0040
CORTEX	Internal Hex Classix	3,3/3,8/4,2/5/6	3,5	DAS_I_0040
CORTEX	Internal Hex Saturn	3,3/3,8/4,2/5/6	3,5	DAS_I_0040
CORTEX	Multi Unit	Universal	Universal	DAS_I_0025
CONEXÃO SISTEMA DE PRÓTESE	Flash	3,5/4,3/5	Universal	DAS_I_0021
CONEXÃO SISTEMA DE PRÓTESE	Torq	3,5/3,75/4	Universal	DAS_I_0021
CONEXÃO SISTEMA DE PRÓTESE	Expand	3,75/4/5	Universal	DAS_I_0021
DENTAURUM	Tiologic	3,3	Small	DAS_I_0130
DENTAURUM	Tiologic	4,8/5,5	Large	DAS_I_0132
DENTEGRIS	SLS-Straight	4.5	4.5	DAS_I_0041
DENTEGRIS	Sinus-Lift	4.5	4.5	DAS_I_0041
DENTIUM	SuperLineImplantium	3,4	3,6	DAS_I_0030
DENTIUM	SuperLineImplantium	3,8	4	DAS_I_0030
DENTIUM	SuperLineImplantium	4,3	4,5	DAS_I_0030
DENTIUM	SuperLineImplantium	4,8	5	DAS_I_0030
DENTIUM	SuperLineImplantium	4,8	6	DAS_I_0030
DENTIUM	SimpleLine II	3,8/4,3	4,8	DAS_I_0037
DENTIUM	SimpleLine II	3,8/4,3	4,8	DAS_I_0074
DENTIUM	SimpleLine II	4,3/4,8	6,5	DAS_I_0096
DENTIS	OneQ-SL	3,9/4,2/4,7/5,2	Regular	DAS_I_0030
DENTIS	OneQ-SL	6/7/8	Wide	DAS_I_0030
DENTIS	s-Clean Tapered	3,7	Mini	DAS_I_0030
DENTIS	s-Clean Tapered	4,1/4,3	Regular	DAS_I_0030

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging



BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
DENTIS	s-Clean Tapered	4,8	Wide	DAS_I_0030
DENTIS	s-Clean Tapered II	3,7	Mini	DAS_I_0030
DENTIS	s-Clean Tapered II	4,1/4,3	Regular	DAS_I_0030
DENTIS	s-Clean Tapered II	4,8	Wide	DAS_I_0030
DENTIS	s-Clean Straight	4,1/4,8	4,1/4,8	DAS_I_0030
DENTIS	s-Clean save	5,5/6	5,5/6	DAS_I_0030
DENTIS	SQ-SL	4/4,5/5	Regular	DAS_I_0030
DENTIS	SQ-SL	6/7/8	Wide	DAS_I_0030
DIO IMPLANTS	External	3,3/3,8	Narrow 3,5	DAS_I_0023
DIO IMPLANTS	External	3,75/4/4,5	Regular 4,1	DAS_I_0024
DIO IMPLANTS	External	5/5,3/5,5/6	Wide 5,1	DAS_I_0061
DIO IMPLANTS	UF	3/3,3	Narrow	DAS_I_0014
DIO IMPLANTS	UF	3,8/4/4,5/5/5,5	qqq	DAS_I_0030
DIO IMPLANTS	Internal OCTA		4.8	DAS_I_0074
DMI Dental Supply	DCI/DSI	3,3/3,5/3,75/4,2/5/6	3,75	DAS_I_0040
DSP BIOMEDICAL	External Hex	3,75/4/5//3,5/3,8/4,3	4,1	DAS_I_0024
DENTAL TECH	Implagic	4,5	4,5 (Blue)	DAS_I_0041
EASY IMPLANT	Mini	3	Mini	DAS_I_0176
EASY IMPLANT	Multi Unit Conical Abutment	Universal	Universal	DAS_I_0025
EASY IMPLANT	Master C	3,5	3,5 (Ocean)	DAS_I_0004
EASY IMPLANT	Master C	4	4 (Ocean)	DAS_I_0004
EASY IMPLANT	Master C	4,5	4,5 (Lilas)	DAS_I_0030
EASY IMPLANT	Master C	5	5 (Lilas)	DAS_I_0030
EASY IMPLANT	Master S	3,3	3,3 (Ocean)	DAS_I_0004
EASY IMPLANT	Master S	3,75	3,75 (Lilas)	DAS_I_0030
EASY IMPLANT	Master S	4,25	4,25 (Lilas)	DAS_I_0030
EASY IMPLANT	Master S	4,75	4,75 (Lilas)	DAS_I_0030
EASY IMPLANT	Master L	3,3	3,3 (Lilas)	DAS_I_0030
EASY IMPLANT	Master L	3,75	3,75 (Lilas)	DAS_I_0030
EASY IMPLANT	Master L	4,25	4,25 (Lilas)	DAS_I_0030
EASY IMPLANT	Master L	4,75	4,75 (Lilas)	DAS_I_0030
ECKERMAN	Hexagon	3/3,5/4/4,5/5	4,1	DAS_I_0024
ELITE MEDICAL	External Connection	3,75	Narrow Platform	DAS_I_0023

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
ELITE MEDICAL	External Connection	4	Regular Platform	DAS_I_0024
ELITE MEDICAL	Conexión Externa	5	Wide	DAS_I_0061
EUROTEKNIKA	Aesthetica	4,1	4,8	DAS_I_0037
EUROTEKNIKA	Aesthetica	4,1	4,8	DAS_I_0074
EUROTEKNIKA	Aesthetica	4,8	6,5	DAS_I_0096
EUROTEKNIKA	Naturactis	3,5	3,4	DAS_I_0004
EUROTEKNIKA	Naturactis	4	3,8	DAS_I_0004
EUROTEKNIKA	Naturactis	4,5	4,3	DAS_I_0004
EUROTEKNIKA	Naturactis	5	4,8	DAS_I_0004
EUROTEKNIKA	Natea	3,6/4,1/4,8	Narrow	DAS_I_0004
EUROTEKNIKA	Natea	3,6/4,1/4,8	Regular	DAS_I_0004
EUROTEKNIKA	Natea	6	Wide	DAS_I_0004
EUROTEKNIKA	Naturall	3,5	Narrow	DAS_I_0004
EUROTEKNIKA	Naturall	4/4,5	Regular	DAS_I_0004
EUROTEKNIKA	Naturall	5	Wide	DAS_I_0005
EUROTEKNIKA	Uneva	3,6/4,1/4,8/6	4,1	DAS_I_0024
EUROTEKNIKA	Multi Unit Tetra	Universal	Universal	DAS_I_0025
GALIMPLANT	Conexión interna	3,5	3,5	DAS_I_0004
GALIMPLANT	Conexión interna	4	4	DAS_I_0004
GALIMPLANT	Conexión interna	5	5	DAS_I_0004
GALIMPLANT	Conexión externa	3,5/4	4	DAS_I_0024
GALIMPLANT	Multi Unit Multiposicion recto	Universal	Universal	DAS_I_0025
GALIMPLANT	Multi Unit Multiposicion angulado	Universal	Universal	DAS_I_0025
GLOBAL D (TEKKA)	InKone Primo	3,5/4/4,5/5	Universal	DAS_I_0152
GLOBAL D (TEKKA)	InKone Universal	3,5/4/4,5/5	Universal	DAS_I_0152
GMI	Frontier	3,3/3,75/4,25	RP 3,3	DAS_I_0040
GMI	Frontier	4,75/5,75	WP 4,3	DAS_I_0041
GMI	Phoenix	3,3/3,75/4	Standard 4,1	DAS_I_0024
GMI	Phoenix	5	Wide 5,1	DAS_I_0061
GMI	Pilar multiestético angulado	Universal	RP 4,8	DAS_I_0025
GT MEDICAL	Best Fit Internal Hex	3,7/4,1/4,3/4,8	Wide	DAS_I_0005
GT MEDICAL	Best Fit External Hex	3,5	Narrow	DAS_I_0023
GT MEDICAL	Best Fit External Hex	4,1	Regular	DAS_I_0024

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging



BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
GT MEDICAL	Best Fit Hexágono Externo	5,1	Wide	DAS_I_0061
GT MEDICAL	Best Fit Internal Octagon	3,7/4,3/4,8	Regular	DAS_I_0037
GT MEDICAL	Best Fit Internal Octagon	3,7/4,3/4,8	Regular	DAS_I_0074
HANH IMPLANT	Tapered Implant	3,5/4,3	3,5 (Purple)	DAS_I_0021
HANH IMPLANT	Tapered Implant	5	5 (Blue)	DAS_I_0022
HANH IMPLANT	Tapered Implant	7	7 (Green)	DAS_I_0124
HANH IMPLANT	MultiUnit Abutment system	Universal	Universal	DAS_I_0025
HIOSSSEN	ETII SA / ETIII SA	3,5	Mini	DAS_I_0029
HIOSSSEN	ETII SA / ETIII SA	4/4,5/5	Regular	DAS_I_0030
HIOSSSEN	ETIII BA	3,5	Mini	DAS_I_0029
HIOSSSEN	ETIII BA	4/4,5/5	Regular	DAS_I_0030
HITEC	Tapered Self Thread	3,3/3,75	3,5	DAS_I_0040
HITEC	Tapered Self Thread	4,2/5	4,5	DAS_I_0041
HITEC	Logic Plus	3,5	3,7	DAS_I_0040
HITEC	Logic Plus	4,3	3,9	DAS_I_0040
IBS	Magic FC	4/4,5/5/5,5/6/6,5	3,8	DAS_I_0030
IBS	N.R. Fix	3/3,5	3,8	DAS_I_0030
IDO IMPLANTS	I Do Implants	3,8/4/4,5/5/5,5/6/7	Universal	DAS_I_0030
IMPLANT DIRECT	Interactive	3,2/3,7	3	DAS_I_0021
IMPLANT DIRECT	Interactive	4,3/5	3,4	DAS_I_0022
IMPLANT DIRECT	Legacy	3,7/4,2	3,5	DAS_I_0040
IMPLANT DIRECT	Legacy	4,7	4,5	DAS_I_0041
IMPLANT DIRECT	Replus/Replant/Reactive	3,5/3,7/4,2	3,5	DAS_I_0026
IMPLANT DIRECT	Replus/Replant/Reactive	4,3/4,7	4,3	DAS_I_0027
IMPLANT DIRECT	Replus/Replant/Reactive	5/5,7	5	DAS_I_0028
IMPLANT DIRECT	SwishActive	3,3	3	DAS_I_0021
IMPLANT DIRECT	SwishActive	4,1/4,8	3,4	DAS_I_0022
IMPLANT DIRECT	Swishplant/ Swishplus	4,1/4,8	4,8	DAS_I_0037
IMPLANT DIRECT	Swishplant/ Swishplus	4,1/4,8	4,8	DAS_I_0074
IMPLANT DIRECT	Swishplant/ Swishplus	4,8/5,7	6,5	DAS_I_0096
IMPLANT GENESIS	AKTIVE	3,5/3,75/4,2/5	Standard	DAS_I_0040
IHDE DENTAL (IMBIODENT)	Bone Level Plus	3,3	3,3	DAS_I_0033
IHDE DENTAL (IMBIODENT)	Bone Level Plus	4,1	4,1	DAS_I_0035

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
IHDE DENTAL (IMBIODENT)	Bone Level Plus	4,8	4,8	DAS_I_0035
INTRA LOCK	IntraHex	3,75/4	3,5	DAS_I_0040
INTRA LOCK	IntraHex	4,75	4,5	DAS_I_0041
INTRA LOCK	Unihex	4	Regular	DAS_I_0024
INTRA LOCK	Unihex	4,75	Wide	DAS_I_0024
JDENTALCARE	JDEvolution/JDEvolution Plus	3,7	3,7	DAS_I_0040
JDENTALCARE	JDEvolution/JDEvolution Plus	4,3/5	4	DAS_I_0040
JDENTALCARE	JDEvolution/JDEvolution Plus	6	5	DAS_I_0040
JDENTALCARE	JD ICON	3,9	3,9	DAS_I_0022
JDENTALCARE	JD ICON	4,3	4	DAS_I_0022
JDENTALCARE	JD ICON	5	4,7	DAS_I_0022
JDENTALCARE	JD ICON PLUS	3,7	3,7	DAS_I_0022
JDENTALCARE	JD ICON PLUS	4,3	4	DAS_I_0022
JDENTALCARE	JD ICON PLUS	5	4,8	DAS_I_0022
KEYSTONE	Internal TiLobe PrimaConnex	3,3/3,5	3,5	DAS_I_0044
KEYSTONE	Internal TiLobe PrimaConnex	4/4,1	4,1	DAS_I_0045
KEYSTONE	Internal TiLobe PrimaConnex	5	5	DAS_I_0046
KEYSTONE	Internal TiLobe Prima Plus	3,5	3,5	DAS_I_0044
KEYSTONE	Internal TiLobe Prima Plus	4,1	4,1	DAS_I_0045
KEYSTONE	Internal TiLobe Prima Plus	5/6	5	DAS_I_0046
KEYSTONE	Restore	3,75/4	RD 4,1	DAS_I_0024
KLOCKNER	Essential Cone	3,5/4/4,5	4,5	DAS_I_0054
KLOCKNER	KL	3,3	Narrow	DAS_I_0023
KLOCKNER	KL	3,7/4,2	Regular	DAS_I_0024
KLOCKNER	KL	4,7	Wide	DAS_I_0061
KLOCKNER	Vega	3,5	NV	DAS_I_0082
KLOCKNER	Vega	4/4,5	RV	DAS_I_0083
KLOCKNER	Vega +	3,6	NV	DAS_I_0082
KLOCKNER	Vega +	4,1/4,6	RV	DAS_I_0083
LASAK	Bioniq	2,9	QN (Amarillo)	DAS_I_0166
LASAK	Bioniq	3,5/4/5	QR(Azul)	DAS_I_0167
LASAK	Multi Unit	Universal	Universal	DAS_I_0168
LEADER	Tixos External Hex	3,3/3,75	4,1	DAS_I_0024

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging



BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
LEADER	Tixos Internal Hex	3,3	3,5	DAS_I_0040
LEADER	Tixos Internal Hex	3,75	4	DAS_I_0040
MEDENTIS (ICX)	Premium/Active Master	3,75	3,75	DAS_I_0125
MEDENTIS (ICX)	Premium/Active Master	4,1	4,1	DAS_I_0125
MEDENTIS (ICX)	Premium/Active Master	4,8	4,8	DAS_I_0125
MEDENTIS (ICX)	Premium/Active Master	3,3	XNP 3,3	DAS_I_0249
MEGAGEN	AnyOne External	3,5	Small 3,5	DAS_I_0023
MEGAGEN	AnyOne External	4	Regular 4,1	DAS_I_0024
MEGAGEN	AnyOne External	4,5	Regular 4,5	DAS_I_0024
MEGAGEN	AnyOne Internal	3,5/4/4,5/5/6/7	General	DAS_I_0030
MEGAGEN	AnyRidge	3,5	Small	DAS_I_0015
MEGAGEN	AnyRidge	4/4,5	Regular	DAS_I_0015
MEGAGEN	AnyRidge	5/5,5	Wide	DAS_I_0015
MEGAGEN	Mini	3/3,4	Mini	DAS_I_0014
MEGAGEN	Cone Abutment	Universal	4,8	DAS_I_0074
MEGAGEN	Multi Unit N Type	Universal	Universal	DAS_I_0025
MICRODENT	Ektos	3,7/4,2	3,5	DAS_I_0040
MICRODENT	System	2,8/3,25	3,5	DAS_I_0003
MICRODENT	Universal	2,8/3,25	3,5	DAS_I_0003
MICRODENT	Universal	3,3/3,5/3,75/4	4,1	DAS_I_0024
MIS	C1	3,3	Narrow	DAS_I_0016
MIS	C1	3,75/4,2	Standard	DAS_I_0017
MIS	C1	5	Wide	DAS_I_0018
MIS	V3	3,9/4,3/5	Standard	DAS_I_0017
MIS	Lance	3,75/4,2	Standard	DAS_I_0024
MIS	M4	3,3	Narrow	DAS_I_0019
MIS	M4	3,75/4,2	Standard	DAS_I_0040
MIS	M4	5/6	Wide	DAS_I_0041
MIS	MultiUnit	Universal	General	DAS_I_0020
MIS	Seven	3,3	Narrow	DAS_I_0019
MIS	Seven	3,75/4,2	Standard	DAS_I_0040
MIS	Seven	5/6	Wide	DAS_I_0041
MOZOGRAU (Ticare)	MG Inhex	3,3	2,3 Mini	DAS_I_0109

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
MOZOGRAU (Ticare)	MG Inhex	3,75/4,25	2,8 Standard	DAS_I_0004
MOZOGRAU	MG Inhex	5	3,8 Maxi	DAS_I_0005
MOZOGRAU	MG Osseous	3,3	3,4 Mini	DAS_I_0003
MOZOGRAU	MG Osseous	3,4/3,75/4,25	4,1 Standard	DAS_I_0024
MOZOGRAU	MG Osseous	5	5 Maxi	DAS_I_0061
MPI	Privilege CM	3,5/4	Regular	DAS_I_0004
MPI	Privilege CM	5	Wide	DAS_I_0005
MPI	Excellence CM	3,5/4	Regular	DAS_I_0004
MPI	Excellence CM	5	Wide	DAS_I_0005
MPI	External Connection HE	3,3/4	4,1	DAS_I_0024
NEOBIOTECH	EB External System	3,5	Narrow	DAS_I_0023
NEOBIOTECH	IS System	4	Regular 4	DAS_I_0030
NEOBIOTECH	IS System	4,5	Regular 4,5	DAS_I_0030
NEOBIOTECH	IS System	5	Wide 5	DAS_I_0030
NEODENT	Cono Morse Mini Pilar	Universal	Regular	DAS_I_0025
NEODENT	Smart HE	3,75/4	4,1	DAS_I_0024
NEODENT	GM	3,5/3,75/4/4,3/5/6	Regular	DAS_I_0186
NEOSS	ProActive Straight/Tapered/Edge	3,5/4/4,5	Standard	DAS_I_0047
NEOSS	ProActive Straight/Tapered/Edge	5/5,5	Standard	DAS_I_0048
NEOSS	ProActive Wide	6	Standard	DAS_I_0048
NEOSS	ProActive Sinus	6,5	Standard	DAS_I_0048
NEOSS	Short Implant	3,5/4/4,5	Standard	DAS_I_0047
NEOSS	Short Implant	5/5,5/6/6,5	Standard	DAS_I_0048
NOBEL BIOCARE	Active	3	Mini 3.0	DAS_I_0159
NOBEL BIOCARE	Active	3,5	Narrow	DAS_I_0021
NOBEL BIOCARE	Active	4,3/5	Regular	DAS_I_0022
NOBEL BIOCARE	Active	5,5	Wide	DAS_I_0124
NOBEL BIOCARE	Branemark	3,3	Narrow	DAS_I_0023
NOBEL BIOCARE	Branemark	3,75/4	Regular	DAS_I_0024
NOBEL BIOCARE	Branemark	5/6	Wide	DAS_I_0061
NOBEL BIOCARE	MultiUnit	4,8	Regular	DAS_I_0025
NOBEL BIOCARE	Replace	3,5	Narrow	DAS_I_0026
NOBEL BIOCARE	Replace	4,3	Regular	DAS_I_0027

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging



BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
NOBEL BIOCARE	Replace	5	Wide	DAS_I_0028
NOBEL BIOCARE	Replace	6	Platform 6	DAS_I_0129
NOBEL BIOCARE	NobelSpeedy	3,3	Narrow	DAS_I_0023
NOBEL BIOCARE	NobelSpeedy	4/5	Regular	DAS_I_0024
NOBEL BIOCARE	NobelSpeedy	5/6	Wide	DAS_I_0061
NORIS MEDICAL	TUFF	3,3/3,75/4,2/5/6	3,75	DAS_I_0040
NORIS MEDICAL	TUFF TT	3,3/3,75/4,2/5/6	3,75	DAS_I_0040
NORIS MEDICAL	ONIX	3,3/3,75/4,2/5/6	3,75	DAS_I_0040
NORIS MEDICAL	CORTICAL	4/5/6.	3,75	DAS_I_0040
NORIS MEDICAL	PTERYCORE	4,2	3,75	DAS_I_0040
NORIS MEDICAL	PTERYFIT	4,2	3,75	DAS_I_0040
NORMON	HI	3,75/4,25/4,75	3,5	DAS_I_0040
NORMON	HE	3,25/3,75/4,25/4,75	4,1	DAS_I_0024
NOVA IMPLANTS	PSI/PCI	3,3/3,75/4,2/5/6	3,75	DAS_I_0040
OSSEOLIFE (ZIACOM)	OEX	3,75/4,25	RP 4,1	DAS_I_0024
OSSTEM	TS	3,5	Mini	DAS_I_0029
OSSTEM	TS	4/4,5/5/6/7	Regular	DAS_I_0030
OSSTEM	US	3,3/3,5	Mini 3,5	DAS_I_0023
OSSTEM	US	3,75/4/4,5	Regular 4,1	DAS_I_0024
OSSTEM	US	5/5,5	Wide 5,1	DAS_I_0061
OSTEOPLUS	SHE	3,45	3,45	DAS_I_0009
OSTEOPLUS	SHE	3,75 / 4	4	DAS_I_0024
OSTEOPLUS	SHI	3,3 / 3,75 / 4,2	3,5	DAS_I_0040
OXY	K1 Line	3,5/4/4,5/5	Regular	DAS_I_0015
OXY	PSK Line	3,5/4/4,5/5	Regular	DAS_I_0015
OXY	MD Line KONE	4,75/4,25/5	Regular	DAS_I_0015
OXY	MD Line Ext	3,75/4,25	Standard	DAS_I_0024
PALTOP	ADVANCED CLASSIC	3,75/4,2/5	Standard	DAS_I_0040
PALTOP	ADVANCED +	3,75/4,2/5	Standard	DAS_I_0040
PALTOP	DYNAMIC	3,75/4,2/5	Standard	DAS_I_0040
PALTOP	DIVA	3,75/4,2/5	Standard	DAS_I_0040
PALTOP	ACTIVE	3,75/4,2/5	Standard	DAS_I_0040
PALTOP	CONICAL ACTIVE	3,75/4,2/5	Standard	DAS_I_0029

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging

BRAND	MODEL	Ø IMPLANT	PLAT-FORM	LIBRARY CODE
PALTOP	Universal Multi-Unit		Universal	DAS_I_0181
PHIBO	BTN Serie 3	3,6	4	DAS_I_0024
PHIBO	BTN Serie 4	4,2	4	DAS_I_0024
PHIBO	TSH Serie 3	3,6	4	DAS_I_0024
PHIBO	TSH Serie 4	4,2	4	DAS_I_0024
PROCLINIC	Cilindrico Externo/Cónico Externo	3,75/4,25//3,5/4	4,1 Estandar	DAS_I_0024
PROCLINIC	Cilindrico Interno/Cónico Interno	3,3/3,75/4,25/5//3,5/4/5	3,5	DAS_I_0040
PROCLINIC	SP Octa	3,3/4,1/4,8	4,8	DAS_I_0037
PROCLINIC	SP Octa	3,3/4,1/4,8	4,8	DAS_I_0074
PROCLINIC	SP Octa	4,8	6,5	DAS_I_0096
PROCLINIC	Aqua CM	3,5/4/5	2,82	DAS_I_0004
RADHEX	PHE	3,5	3,5	DAS_I_0023
RADHEX	PHE	4/4,5/5	4,1	DAS_I_0024
RADHEX	PHI	3,75	3,5	DAS_I_0040
RADHEX	PHI	4,5/5	4,5	DAS_I_0041
SEWON MEDIX	IH2 SLA SYSTEM	3,5	Mini	DAS_I_0029
SEWON MEDIX	IH2 RBM SYSTEM	3,5	Mini	DAS_I_0029
SEWON MEDIX	IH2 SLA SYSTEM	3,5/4/4,5/5	Regular	DAS_I_0030
SEWON MEDIX	IH2 RBM SYSTEM	3,5/4/4,5/5	Regular	DAS_I_0030
SEWON MEDIX	IH SYSTEM	3,5/4/4,5/5	Universal	DAS_I_0025
SGS	P1	3,2/3,75/4,2/5/6	3,75	DAS_I_0040
SGS	P7	3,2/3,75/4,2/4,5/5/6	3,75	DAS_I_0040
SIGNO VINCES	Duo	4,6	4,1	DAS_I_0024
SIGNO VINCES	Inttegra	3,75/4	4,1	DAS_I_0024
SIGNO VINCES	Compact	4,5	CM3,8	DAS_I_0004
SIGNO VINCES	Duocon	3,8	CM3,8	DAS_I_0004
SIGNO VINCES	Duocon	4,6/5,5	CM4,6	DAS_I_0005
SIGNO VINCES	Infra	3,3/3,8/4,6	CM	DAS_I_0004
SIN IMPLANTS	HI SW	3,8	3,8	DAS_I_0039
SIN IMPLANTS	Tryon	3,25/3,75/4	4,1	DAS_I_0024
SIN IMPLANTS	Tryon CO	4	4,1	DAS_I_0024
SIN IMPLANTS	Revolution	3,25/3,75/4	4,1	DAS_I_0024
SIN IMPLANTS	Stylus	4	4,1	DAS_I_0024

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging



BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
SOUTHERN IMPLANTS	Tri-Nex	3,5	3,5	DAS_I_0026
SOUTHERN IMPLANTS	Tri-Nex	4,3	4,3	DAS_I_0027
SOUTHERN IMPLANTS	Tri-Nex	5	5	DAS_I_0028
SOUTHERN IMPLANTS	Tri-Nex	6	6	DAS_I_0129
SOUTHERN IMPLANTS	IT Connection	3,3/4/4,1/4,9/5	4,8	DAS_I_0037
SOUTHERN IMPLANTS	IT Connection	4,9/5/6	6,5	DAS_I_0096
SOUTHERN IMPLANTS	External Hex	3,25	3,4	DAS_I_0003
SOUTHERN IMPLANTS	External Hex	3,75/4	4,1	DAS_I_0024
SOUTHERN IMPLANTS	Internal Hex	3,75/4,2/5	3.6	DAS_I_0040
SOUTHERN IMPLANTS	Compact Conical MultiUnit	4,8	4,8	DAS_I_0025
SOUTHERN IMPLANTS	Deep Conical	3	2,45	DAS_I_0109
SOUTHERN IMPLANTS	Deep Conical	3,5/4	2,95/3,1	DAS_I_0004
SOUTHERN IMPLANTS	Deep Conical	5	4,1	DAS_I_0005
SOUTHERN IMPLANTS	Provata	4/5/6	3.6	DAS_I_0040
STERIOSS	HexLoc	3,25	3,3	DAS_I_0023
STERNGOLD	Stern Ex	3,75/4/5	4,1	DAS_I_0024
SYBRON	Endopore	4,1	4,1	DAS_I_0024
SWEDEN MARTINA	Premium khono	3,3	3,3	DAS_I_0031
SWEDEN MARTINA	Premium khono	3,8	3,8	DAS_I_0032
SWEDEN MARTINA	Premium khono	4,25	4,25	DAS_I_0065
SWEDEN MARTINA	Prama	3,8/4,25/5	3,4	DAS_I_0031
SWEDEN MARTINA	Shelta	3,8	3,8	DAS_I_0032
SWEDEN MARTINA	Shelta	4,25	4,25	DAS_I_0032
SWEDEN MARTINA	Shelta	5	5	DAS_I_0032
STRAUMANN	Standard Plus	3,3	NNC 3,5	DAS_I_0160
STRAUMANN	Bone Level	3,3	NC 3,3	DAS_I_0033
STRAUMANN	Bone Level	4,1/4,8	RC4,1/4,8	DAS_I_0035
STRAUMANN	BLX	3,5/3,75/4/4,5	RB (Regular Base)	DAS_I_0207
STRAUMANN	BLX	5/5,5/6,5	WB (Wide Base)	DAS_I_0208
STRAUMANN	Tissue Level	3,3/4,1/4,8	Regular 4,8	DAS_I_0037
STRAUMANN	Tissue Level	4,8	Wide 6,5	DAS_I_0096
STRAUMANN	Synocta	Universal	Regular 4,8	DAS_I_0074
STRAUMANN	Screw Retained	Universal	NC/RC	DAS_I_0101

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging

BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
STRAUMANN	Tissue Level	3,3	3,5	DAS_I_0160
TBR	Hexágono Externo	3,5	Narrow	DAS_I_0023
TBR	Hexágono Externo	5	Wide	DAS_I_0058
TITANIUMFIX	Bfix	3,5/4	Regular	DAS_I_0004
TITANIUMFIX	Bfix	4,5/5	Larga	DAS_I_0005
TREE-OSS	Simple	3,3/3,75/5	3,75	DAS_I_0040
TREE-OSS	Hexagono Externo	3,3	3,5 Rosa	DAS_I_0023
TREE-OSS	Hexagono Externo	3,75/4	4,1 Amarillo	DAS_I_0024
TREE-OSS	Hexagono Externo	5	5,1 Azul	DAS_I_0061
TREE-OSS	Tri-Channel Interno	3,5	3,5 Rosa	DAS_I_0026
TREE-OSS	Tri-Channel Interno	4,3	4,3 Amarillo	DAS_I_0027
TREE-OSS	Tri-Channel Interno	5	5,0 Azul	DAS_I_0028
TREE-OSS	MULTI UNIT	Universal	Universal	DAS_I_0025
TRI DENTAL IMPLANTS	TriVent	3,75/4,1/4,7	3,5	DAS_I_0040
TRINON	Q2	3,5/3,75/4,5	4	DAS_I_0024
TRINON	QK	4	4,8	DAS_I_0037
TRINON	QK	4	4,8	DAS_I_0074
UFIT	Gt2	3,5	Mini	DAS_I_0004
UFIT	Gt2	4/4,5	Regular	DAS_I_0005
UFIT	Gt2	5	Wide	DAS_I_0005
UFIT	Gt2	5,5/6/6,5/7	Ultrawide	DAS_I_0005
UFIT	Gt2	3,5	Mini	DAS_I_0004
UFIT	Nt2	3,5	Mini	DAS_I_0004
UFIT	Nt2	4/4,5	Regular	DAS_I_0005
UFIT	Nt2	5	Wide	DAS_I_0005
UFIT	Nt2	5,5/6/6,5/7	Ultrawide	DAS_I_0005
VULKAN IMPLANTS	INHex		3,75	DAS_I_0040
XIVE	S	3,4	3,4	DAS_I_0038
XIVE	S	3,8	3,8	DAS_I_0039
XIVE	S	4,5	4,5	DAS_I_0085
XIVE	XIVE	5,5	5,5	DAS_I_0086
YES IMPLANT	SSYSTEM	3,3/3,5	Narrow	DAS_I_0030
YES IMPLANT	SSYSTEM	4/4,5	Regular	DAS_I_0030

NOTE:

DAS = Dynamic Abutment Solutions

I = Intraoral/Extraoral

8/10/12 = Scanbody heights

0001/0002/0003/... = Compatibility

Premium = Premium will only appear if you have purchased the DAS Premium Library.

GH= Gingival Height

CH= Cement Height

R= Rotational/Non-Engaging

NR= Non Rotational/Engaging



BRAND	MODEL	Ø IMPLANT	PLATFORM	LIBRARY CODE
YES IMPLANT	SSYSTEM	5/5,5	Wide	DAS_I_0030
YES IMPLANT	SSYSTEM	6/6,5/7/8	BWide	DAS_I_0030
ZIMMER	SwissPlus	3,7/4,1/4,8	4,8	DAS_I_0037
ZIMMER	SwissPlus	3,7/4,1/4,8	4,8	DAS_I_0074
ZIMMER	Zimmer Type (Internal Hexagon)	3,7/4,1	3,5	DAS_I_0040
ZIMMER	Zimmer Type (Internal Hexagon)	4,7	4,5	DAS_I_0041
ZIMMER	ScrewVent	3,7/4,1	3,5	DAS_I_0040
ZIMMER	ScrewVent	4,7	4,5	DAS_I_0041
ZIMMER	ScrewVent	6	5,7	DAS_I_0080
ZIMMER	Tapered	Universal	Universal	DAS_I_0205